Retraction System for Minimally Invasive CABG and Transapical Valve Implantation



Aesculap Surgical Technologies - Surgical Instruments



Jochen Cremer, M.D. Department of Cardiovascular Surgery, University Hospital Kiel, Germany

HSR Proceedings in Intensive Care and Cardiovascular Anesthesia 2011; 3(4): 249-253

The MIDCAB approach in its various dimensions

J. Cremer, J. Schoettler, A. Thiem, C. Grothusen, G. Hoffmann Department of Cardiovascular Surgery, Christian-Albrechts-University of Kiel, Germany

ABSTRACT

Introduction: Minimally invasive direct coronary artery bypass grafting (MIDCAB) offers arterial revascularization of the left anterior descending coronary artery especially in lesions unsuitable for percutaneous coronary interventions. By avoidance of sternotomy and cardiopulmonary bypass its invasivity is less than that of conventional bypass surgery.

Methods: A literature search of all published minimal invasive direct coronary artery bypass grafting studies was performed for the period from January 1995 through November 2011. Additionally, the authors reviewed their experience in more than one thousand patients treated by minimal invasive direct coronary artery bypass grafting within the last 14 years at their institution.

Results: Early mortality ranged from 1.2 to 1.3 %. Midterm mortality ranged up to 3.2 %. At 6-month follow up 3.6 % of grafts were occluded and 7.2 % had a significant stenosis which resulted in target vessel revascularization in 3.3 % of cases. The conversion rate to sternotomy or cardiopulmonary bypass ranged between 1.2 and 6.2 %.

Conclusions: In the past MIDCAB was predominantly used in patients with isolated lesions of the left anterior descending coronary artery. In combination with percutaneous interventions it provides an attractive option for full revascularization in multi vessel disease especially in older patients with significant comorbidities. Overall minimal invasive direct coronary artery bypass grafting is associated with few perioperative complications and with high graft patency rates in the mid-term and long-term course.

Keywords: coronary artery disease, arterial revascularization, MIDCAB, hybrid procedure, minimally invasive cardiac surgery.

Presented at the Roland Hetzer International Cardiothoracic and Vascular Surgery Society 1st Expert Forum

INTRODUCTION

After the first description of a beating heart anastomosis between the mammary artery and a coronary artery by Kollesov (1) this concept wasn't established in clinical reality for decades. Since its reintroduction in the 1990s by Calafiore (2) and Subramanian (3) minimally invasive direct coronary artery bypass grafting (MIDCAB) into the spectrum of surgical revascularization this technique has been further developed. Today it is a representative part of the cardiac surgical program in some institutions.

The procedure is regarded as less invasive than conventional coronary artery bypass (CABG) or off pump revascularization (OPCAB) and ranges in its invasivity close to percutaneous coronary interventions (PCI). MIDCAB can be regarded as an

Jochen Cremer, M.D. Department of Cardiovascular Surgery, University Hospital Kiel, Germany

alternative treatment option besides the well established treatment procedures of revascularization.

MIDCAB has been mainly used in patients with proximal stenosis of the left anterior descending artery (LAD). In these patients interventional treatment with PCI appears risky or impossible due to complex lesions, close relationship to the main stem or other coronary arteries or total occlusion of the target vessel. In other patients repeated interventions at the LAD remained without longstanding success.

Besides its original intention for revascularization of the LAD, MIDCAB can be a useful part of hybrid procedures in patients with multi vessel disease where a major coronary surgery procedure would not be well tolerated. Although multi vessel disease is a predictor of elevated mortality after CABG, the MIDCAB procedure can be performed as stand alone procedure with acceptable results in midterm morbidity and mortality, even though formally incomplete revascularization may remain (4). Several studies and our own experience proved that in selected patients with main stem stenosis or multi vessel disease MIDCAB can be safely performed (5). Complete revascularization can be achieved by a hybrid approach with accompanying PCI (5-6). Interestingly, we found in our own collective that a significant part of the planned PCI has not been performed since the MIDCAB resulted in ongoing freedom from angina. This is particularly true in multi vessel disease where the LAD is the dominating vessel.

As a rare indication we treated two children aged twelve and thirteen who had undergone a switch operation as babies for transposition of the greater arteries. Later they developed significant stenoses of the LAD and were successfully treated by MIDCAB with an uneventful postoperative course and restored coronary circulation.

Recently, we demonstrated that MIDCAB can be performed in octogenarians with satisfying midterm results. Especially these patients benefit from the marginal surgical trauma, less anesthesia and short ventilation duration and intensive care unit stay (7). A small subgroup of patients received MIDCAB who suffered from severe bleeding complications due to dual antiplatelet medication after PCI with drug eluting stents.

METHODS

The patient is placed under general anesthesia and is intubated with a doublelumen tube. Standard monitoring includes 5-lead electrocardiogram (ECG), arterial blood pressure, central venous pressure, nasopharyngeal temperature and peripheral oxygen saturation. The main features of the procedure consist of a small left anterior thoracotomy in the fourth or fifth intercostal space. The pleural cavity is opened and to facilitate harvesting of the graft the left lung is deflated. The left internal thoracic artery (LITA) is harvested as a pedicle up to the level of the first rib. A thoracic wall tilting device attached to the upper rib is useful to reach the proximal parts of the LITA. After completion of harvest the LITA is dissected distally at the level of the sixth intercostal space. Graft perfusion is controlled and local application of papaverine as vasodilating agent can optimize graft flow.

After opening the pericardium and placement of fraction sutures the target region is partially immobilized using a special mechanical stabilizer. Before occlusion of the LAD 100 IU/kg of unfractioned heparin are administered. In almost every case preconditioning is used by temporary occlusion of the LAD with encircling 4-0 tourniquets sutures over silicone tubing. After a short period of reperfusion the target vessel is incised. The field is kept free of blood by a humidified blower using fluid and carbon dioxide. The anastomosis between LITA and LAD is performed under direct vision using an 8/0 monofilament suture, whereas the first five stitches are performed in distance. In contrast to other groups videoscopy as well as intracoronary shunts are not applied.

Following completion of the anastomosis the pedicle is fixed with two sutures to avoid torsion of the graft. After restoration of blood flow heparin is antagonized with protamine. The initially mobilized prepericardial fat is fixed to the medial pericardium to cover the distal past of the IMA course. An infero-posterior chest tube is placed in the left pleural cavity. Closing of the thoracotomy including 2 transcostal sutures has to be performed with meticulous care to avoid contact of the graft with the chest wall or herniation of mediastinal structures. On table extubation and a short intensive care unit stay is intended in all patients.

RESULTS

We started the MIDCAB program in our institution in 1998. Meanwhile over 1000 patients had been operated with a rather low conversion rate of 1.2 %. The majority of patients had a single vessel disease (55.6 %) and half of the patients were pre-treated with PCI or had a complete occlusion of the LAD. The overall 30 day mortality until end of 2010 was 1.2 % (11/914), whereas a substantial amount of these patients was accepted with a high risk profile or for an ultima ratio approach.

Perioperative myocardial infarctions or cerebrovascular accidents were comparably rare events. Rethoracotomy for bleeding complications was necessary in 1.3 %.

DISCUSSION

When comparing PCI and MIDCAB several studies and a meta-analysis showed that overall mortality and myocardial infarction rates are similar, but MIDCAB is associated with significantly lower rates of repeat target vessel revascularizations (8-11). It has to be taken into account that these studies compared MIDCAB with bare metal stents. Recently Thiele and co-workers (12) published a randomized trial with a non-inferiority design which compared PCI of the LAD with drug eluting stents and MIDCAB. The results were comparable concerning death and myocardial infarction but the rate of target vessel revascularization was still higher in the PCI group.

At our institution we apply in almost every case ischemic preconditioning (5 min. of ischemia followed by 2 min. of reperfusion). It's potentially beneficial effect on myocardial performance during and after ischemia is discussed controversially in clinical trials and experimental investigations (13-14). Nevertheless the temporary occlusion anticipates the expectable events during anastomosis like ST-elevation, rhythm disturbances and hemodynamic compromise. These events are more likely in partially occluded target vessels and closeness to other coronary vessels.

Our results come along with published data from other large series of MIDCAB (15). The perioperative mortality compares favourably with the 2.1 % mortality for offpump and 5.9 % mortality for on-pump single bypass grafting as reported in the German Society for Thoracic and Cardiovascular Surgery in 2010 (16). Nevertheless, long-term follow-up results and randomized prospective trials are needed to compare MIDCAB surgery with competing revascularization methods.

Jochen Cremer, M.D. Department of Cardiovascular Surgery, University Hospital Kiel, Germany

A meta-analysis of Kettering in 2008 summarizes the results of seventeen outcome studies of MIDCAB (17). Early and late death rates were 1.3 and 3.2 %. At 6-month follow up 3.6 % grafts were occluded and 7.2 % had a significant stenosis which resulted in 3.3 % in a target vessel revascularization. The incidence of myocardial infarction was below 1 %. With respect to the initial learning curve one could assume that the actual results of MIDCAB change for a better. For example, the perioperative mortality in our institution declined to 0.7 % within the last five years. A recent trial showed a graft patency rate of 100 % after 12 months. These results invalidate the initial concerns about the suboptimal quality of beating heart anastomosis via a limited surgical access. Nevertheless these good results are highly dependent on an experienced surgeon.

Limitations. In a small amount of cases (1.2 to 1.7 % literature and own results) the MIDCAB procedure can't be completed as planned. In some cases the LAD is not to be found due to intramuscular course or excessive epicardial fat tissue or it is situated too laterally to be reached with the LIMA graft. We saw this repeatedly in patients with long standing atrial fibrillation and resulting atrial and right ventricular enlargement. As a result we perform liberally preoperative CT Scans in these patients to judge the reachability in advance. In some cases it was possible to elongate the LITA with a venous interponate. The same strategy of prolongation can be used for injured LITAs for example in local dissection or electrocautery damage.

Malignant rhythm disturbances occur rarely under epicardial manipulation. If these are not of self-limitating nature cardioversion or defibrillation may be neces sary. To prevent these events we advocate keeping the patient on high potassium level and sufficient mean arterial pressure. Additional administration of magnesium may be useful. In high risk patients amiodarone is prophylactically given. Since we repeatedly observed significant bradycardia during operation an epicardial pacing should become quickly available, when necessary. In addition external defibrillation adhesive pads should be used in high risk patients.

In situation with critical ischemia an intra-luminal shunt is placed.

In the early era of MIDCAB, adhesions of the LITA to the thoracic chest wall occurred sporadically. This can be associated with life-threatening bleedings due to disruption of the anastomotic site.

Due to events like this, closing of the operation site in several layers with especially covering the distal IMA course is of paramount importance. By means of that chest wall hernias can also be prevented.

CONCLUSION

Overall we regard MIDCAB revascularization as a very valuable concept within the different options of surgical revascularization. From the beginning it could be proven that the MIDCAB approach is safe and efficient. By avoidance of sternotomy and heart-lung machine it represents a gentle way of effective revascularization of the most important coronary vessel with the best possible graft.

MIDCAB is not only applicable in isolated lesions of the LAD. Hybrid concepts are of growing interest in serving older patients with significant comorbidities.

MIDCAB is associated with a low level of complications perioperatively and with high graft patency rates in the mid-term and long-term course.

References

- Kolessov VI. Mammary artery-coronary artery anastomosis as method of treatment for angina pectoris. J Thorac Cardiovasc Surg. 1967;54:535-44.
- Calafiore AM, Di Giammarco G, Teodori G, et al. Left an terior descending coronary artery grafting via left anterior small thoracotomy without cardiopulmonary bypass. Ann Thorac Surg. 1996;61:1658-65.
- Subramanian VA, Sani G, Benetti FJ, et al. Minimally inva sive coronary artery bypass surgery: a multicenter report of preliminary clinical experience. Circulation. 1995;92:1645.
- Lichtenberg A, Klima U, Paeschke H. Impact of multivessel coronary artery disease on outcome after isolated minimally invasive bypass grafting of the left anterior descending artery. Ann Thorac Surg. 2004;78:487-91.
- Holzhey DM, Jacobs S, Mochalski M, et al. Minimally invasive hybrid coronary artery revascularization. Ann Thorac Surg. 2008;86:1856–60.
- Wittwer T, Cremer J, Boonstra P, et al. Myocardial hybrid revascularization with minimally invasive direct coronary atery bypass grafting combined with coronary angioplasty: preliminary results of a multicenter study. Heart. 2000;83:58–63.
- 7. Hoffmann G, Fraund-Cremer S, Petzina R, et al. Ergebnisse der MIDCAB-Operation bei Patienten über 80 Jahren. J Clin Res Cardiol 100. 2011 April;1498.
- Jaffery Z, Kowalski M, Weaver WD, et al. A meta-analysis of randomized control trials comparing minimally invasive direct coronary artery bypass grafting versus percutaneous coronary intervention for stenosis of the proximal left anterior descending artery. Eur J Cardiothorac Surg. 2007;31:691–7.
- Fraund S, Hermann G, Witzke A, et al. Midterm follow-up after minimally invasive direct coronary artery bypass grafting versus percutaneous coronary intervention techniques. Ann Thorac Surg. 2005; 79:1225-31.
- Diegeler A, Thiele H, Falk V, et al. Comparion of stenting with minimally invasive bypass surgery for stenosis of the left anterior descending artery. N Engl J Med. 2002;347:561-5.
- Cisowski M, Morawski W, Drzewiecki J, et al. Integrated minimally invasive direct coronary artery bypass grafting and angioplasty for coronary artery revascularization. Eur J Cardiothorac Surg 22. 2002:261–5.
- 12. Thiele H, Neumann-Schniedewind P, Jacobs S, et al. Randomized comparison of minimally invasive direct coronary artery bypass surgery versus sirolimus-eluting stenting in isolated proximal left anterior descending coronary artery stenosis. J Am Coll Cariol. 2009;53:2324-31.
- 13. Halkos ME, Kerendi F, Corvera JS. Myocardial protection with postconditioning is not enhanced by ischemic preconditioning. Ann Thorac Surg. 2004;78:961–9.
- 14. Ghosh S, Galianes M. Protection of the human heart with ischemic preconditioning during cardiac surgery: role of cardiopulmonary bypass. J Thorac Cardiovasc Surg. 2003;126: 33-42.
- Holzhey DM, Jacobs S, Mochalski, M, et al. Seven-Year Follow-up after minimally invasive direct coronary artery bypass: experience with more than 1300 patients. Ann Thorac Surg. 2007;83:108– 114.
- Gummert JF, Funkat AK, Beckmann A, et al. Cardiac surgery in germany during 2010: a report on behalf of the german society for thoracic and cardiovascular surgery. Thorac Cardiov Surg. 2011;59: 259–67.
- Kettering K. Minimally invasive direct coronary artery bypass grafting: a meta analysis. J Cardiovasc Surg. 2008;49:793-800.

Cite this article as: Cremer J, Schoettler J, Thiem A, Grothusen C, Hoffmann G. The MIDCAB approach in its various dimensions. HSR Proceedings in Intensive Care and Cardiovascular Anesthesia. 2011;3(4): 249–253

Source of Support: Nil. Conflict of interest: None declared.



As one of the established surgical approaches for myocardial revascularization, MIDCAB via a left anterior minithoracotomy is especially effective for grafting the LAD in situations of one vessel disease, when PCI is regarded as less effective, critical or technically not applicable.

As the surgical part of a hybrid approach, MIDCAB-LAD grafting may also be an attractive alternative for use in combination with PCI techniques for complete revascularization in multi-vessel disease.

Even though, MIDCAB instrumentation has been developed over the last 15 years, major improvements could be achieved using TAVI-MIDCAB Retraction System from Aesculap.



Gentle spreading of the retractor facilitates favorable access to the pericardium with equal distribution of the mechanical stress due to the flexibility of the rotating double blades which conform to the individuals specific rib anatomy.



Selecting different blade lengths supports the concept of individual optimization of the surgical access, depending on the extent of fat tissue and the general thickness of the thoracic wall. Along with pericardial stay sutures, the distance to the anterior heart and the left ventricular apex can be kept very small, through use of the retractor.



This combination will make even a very distal IMA-anastomosis in MIDCAB procedures possible and gives an easily accessible exposure of the LV apex.





The existence of various fixation points combined with the wide range of the mechanical stabilizer, provide optimal immobilization of the target LAD segment allowing free access for the surgeon as well as easy handling for the surgical assistant.

Fixation points





Thus, a precise anastomosis is successfully performed with very little tissue trauma.



Due to the equal distribution of pressure along the opening of the retractor, the skin integrity remains preserved until the end of the procedure.



Even in complex cases with increased surgical risks, the use of the TAVI-MIDCAB Retractor System provides small incisions which can allow for faster primary wound healing.



MIDCAB Retractor complete

FC525R

consisting of: FC520R Retractor FC521R MIDCAB Adapter FC522R MIDCAB Stabilizer FC530R-FC534R various double blades

Features of the newly designed TAVI-MIDCAB retractor system:

- Improved visibility and access to surgical site due to:
 - Low profile design
 - Flexible, interchangeable blades
 - Multitude of stabilizer mounting positions
 - Suction capability
- Reduced tissue and rib trauma due to flexible blades
- Adaptable to most patients through use of interchangeable blades

With these technical features, this retractor is also very applicable for Transapical TAVI (transcatheter aortic valve implantation) procedures, which are generally performed in elderly comorbid patients, whose tissue sensitivity and consistency is already severely impaired. Retractor only



RETRACTOR only, double joint





MIDCAB Adapter



MIDCAB Adapter for FC522R



MIDCAB Stabilizer

FC522R

The double joint retractor allows for individual patient adjustment in a low profile design. Stabilizer can be connected for suction to clear surgical site.

Double Blades (sold individually)



FC530R

Double Blade, 40 mm



FC531R Double Blade, 50 mm



Double Blade, 60 mm



FC533R Double Blade, 70 mm



FC534R

Double Blade, 80 mm

The rotating double blades conform to the changing dimensions of the ribs in the opening process of the retractor. Various blade lengths are available according to the thickness of the chest wall including the option of different lengths to accommodate the breast of female patients.



All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogs and price lists and to take legal actions.

 $^{\odot}2015$ AESCULAP. ALL RIGHTS RESERVED. PRINTED IN THE USA. Aesculap is an equal opportunity employer

Aesculap, Inc. | 3773 Corporate Parkway | Center Valley, PA | 18034 Phone 800-282-9000 | Fax 610-791-6886 | www.aesculapusa.com