

proGAV 2.0

## MRI SAFETY INFORMATION



MR Conditional

Non-clinical testing demonstrated that the proGAV 2.0 valve is MR Conditional. A patient with the proGAV 2.0 valve can be safely scanned immediately after implantation in an MR system meeting the following conditions:

- Static magnetic field of 3-Tesla or less
- Maximum spatial gradient magnetic field of 1,400-Gauss/cm (14.0 T/m) or less
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4-W/kg (First Level Controlled Operating Mode)
- Do not take the *proGAV 2.0* tools into the MR environment. They are MR Unsafe.

Under the scan conditions defined above, the *proGAV 2.0* valve is expected to produce a maximum temperature rise of +3°C after 15 minutes of continuous scanning.

## COMPATIBILITY WITH DIAGNOSTIC PROCEDURES

MRI examinations with field strengths of up to 3.0 tesla and CT examinations can be carried without endangering or impairing the functionality of the shunt. The *proGAV 2.0* is MR Conditional (ASTM F2503-13). The setting of the *proGAV 2.0* will not change when subjected to an MRI of 1.5 T or 3.0 T. All components are visible via X-ray. The provided catheters are MRI Safe. Reservoirs, deflectors and connectors are MR Conditional.

For additional proGAV 2.0 MRI Safety Information, including artifact information, Warnings and Precautions see product IFU [SOP-AIC-5001399](#).