



ELAN[®] 4 Power System

Tool Catalog

ENVISION THE EXTRAORDINARY[™]



AESCULAP[®]

ELAN[®] 4 Power System

Tool Catalog

Table of Contents

› Sagittal Saw Blades (Oscillating Saw Blades)	3-7
› Reciprocating Saw Blades	8-12
› Transversal Saw Blades	13-14
› ELAN 4 Burrs	15-23

Signs and Symbols

 › Usable Working Length (e.g. of Sawblades)	 › Incision Depth	 › Packed Sterile\Unsterile
 › Burr/Blade Length	Ø › Diameter of Burr	② › Tool for Single Use
 › Burr/Blade Width	TC › Tungsten Carbide Burr	» › For Use with Handpiece
→ ← › Burr/Blade Thickness	 › Packing Unit	 › Page
 › Incision Width: Different from Blade Thickness in Sawblades with Off-Set		

All measurements and technical data listed in this brochure are approximate, as they are subject to normal manufacturing tolerances and technical modifications.


Indications for USE:

The ELAN 4 Tools are intended for high speed cutting, sawing and drilling of bone in the fields of Spine, ENT, Neuro and Maxillofacial Surgery.


Sagittal Saw Blades (Oscillating Saw Blades)



GA836



GA831




ELAN[®] 4 Power System


Sagittal Saw Blades (Oscillating Saw Blades)

5 sterile




GP491SU	14 mm	13 mm	0.3 mm	0.3 mm	» GA836 
GP492SU	15 mm	5 mm	0.3 mm	0.3 mm	
GP493SU	20 mm	5 mm	0.3 mm	0.3 mm	
GP494SU	20 mm	10 mm	0.3 mm	0.3 mm	
GP495SU	20 mm	15 mm	0.3 mm	0.3 mm	
GP496SU	25 mm	5 mm	0.3 mm	0.3 mm	
GP497SU	25 mm	12 mm	0.3 mm	0.3 mm	

GP491SU




GP496SU



GP492SU




GP497SU




GP493SU



GP494SU



GP495SU






GP501SU	15 mm	5 mm	0.5 mm	0.5 mm
GP502SU	15 mm	10 mm	0.5 mm	0.5 mm
GP504SU	20 mm	5 mm	0.5 mm	0.5 mm
GP505SU	20 mm	10 mm	0.5 mm	0.5 mm
GP514SU	25 mm	5 mm	0.5 mm	0.7 mm
GP515SU	25 mm	10 mm	0.5 mm	0.7 mm
GP516SU	25 mm	15 mm	0.5 mm	0.7 mm
GP518SU	35 mm	10 mm	0.5 mm	0.5 mm
GP521SU	35 mm	5 mm	0.5 mm	0.5 mm
GP522SU	35 mm	10 mm	0.5 mm	0.7 mm
GP523SU	35 mm	15 mm	0.5 mm	0.7 mm
GP524SU	35 mm	20 mm	0.5 mm	0.7 mm


» GA831



GP501SU




GP502SU




GP504SU




GP505SU



GP514SU



GP515SU




GP516SU




GP518SU



GP521SU




GP522SU



GP523SU



GP524SU



ELAN[®] 4 Power System


Sagittal Saw Blades (Oscillating Saw Blades)

1 sterile




GP535SU	50 mm	10 mm	0.5 mm	0.8 mm
GP536SU	50 mm	15 mm	0.5 mm	0.8 mm
GP537SU	50 mm	20 mm	0.7 mm	0.9 mm
GP538SU	50 mm	25 mm	0.7 mm	0.9 mm
GP539SU	50 mm	30 mm	0.7 mm	0.9 mm

»GA831



GP535SU




GP536SU




GP537SU



GP538SU



GP539SU






GP511SU	35 mm	10 mm	0.5 mm	0.5 mm	10 mm
GP526SU	51 mm	15 mm	0.5 mm	6 mm	20 mm
GP527SU	51 mm	15 mm	0.5 mm	7 mm	20 mm
GP528SU	51 mm	15 mm	0.5 mm	8 mm	20 mm
GP529SU	51 mm	15 mm	0.5 mm	9 mm	20 mm
GP530SU	51 mm	15 mm	0.5 mm	10 mm	20 mm
GP531SU	51 mm	15 mm	0.5 mm	11 mm	20 mm
GP532SU	51 mm	15 mm	0.5 mm	12 mm	20 mm


» GA831



GP511SU



GP526SU -
GP532SU




ELAN® 4 Power System


Reciprocating Saw Blades



GA837



GA832






GP542SU	13 mm	0.3 mm	0.3 mm
GP543SU	20 mm	0.3 mm	0.3 mm
GP544SU	13 mm	0.3 mm	0.3 mm
GP545SU	20 mm	0.3 mm	0.3 mm


» GA837




GP542SU



GP543SU



GP544SU



GP545SU




ELAN[®] 4 Power System

Reciprocating Saw Blades




GP551SU	84 mm	25 mm	0.40 mm	0.60 mm	
GP561SU	104 mm	25 mm	0.40 mm	0.60 mm	
GP552SU	85 mm	27 mm	0.25 mm	0.35 mm	
GP553SU	85 mm	27 mm	0.40 mm	0.60 mm	
GP562SU	105 mm	27 mm	0.25 mm	0.35 mm	
GP563SU	105 mm	27 mm	0.40 mm	0.60 mm	
GP554SU	91 mm	33 mm	0.40 mm	0.60 mm	
GP564SU	111 mm	33 mm	0.40 mm	0.60 mm	

GP551SU



GP561SU



GP552SU

GP553SU




GP562SU

GP563SU



GP554SU



GP564SU





GP565SU

115 mm

37 mm

0.4 mm

0.6 mm

» GA832

GP550SU

73 mm

15 mm

0.4 mm

0.6 mm

GP575SU

79 mm

17 mm

0.4 mm

0.6 mm

4 mm

GP576SU

86 mm

24 mm

0.4 mm

0.6 mm

6 mm



GC615SU

23 mm

0.25 mm

0.35 mm

GC633SU

23 mm

0.25 mm

0.35 mm




GP559SU

70 mm


GP569SU

90 mm


GP565SU




GP550SU




GP575SU



GP576SU




GC615SU



GC633SU



GP559SU



GP569SU




ELAN[®] 4 Power System

Reciprocating Saw Blades




GP557SU	70 mm	10.0 mm	» GA832 
GP558SU	70 mm	12.5 mm	
GP567SU	90 mm	10.0 mm	


GP557SU



GP558SU




GP567SU



Transversal Saw Blades




GA833



ELAN[®] 4 Power System

Transversal Saw Blades




GP581SU	8 mm	4 mm	0.3 mm	0.4 mm	GA833 
GP583SU	10 mm	10 mm	0.5 mm	0.7 mm	
GP584SU	15 mm	12 mm	0.5 mm	0.7 mm	
GP587SU	11 mm	7 mm	0.6 mm	0.7 mm	
GP588SU	11 mm	12 mm	0.6 mm	0.7 mm	
GP591SU	8 mm	25 mm	0.5 mm	0.7 mm	
GP594SU	8 mm	12 mm	0.5 mm	0.7 mm	
GP595SU	14 mm	12 mm	0.5 mm	0.7 mm	
GP598SU	37 mm	25 mm	0.5 mm	0.7 mm	

GP581SU




S Transversal
to be used with
GP580R

GP583SU




M Transversal
to be used with
GP582R & GP586R

GP584SU




M Transversal
to be used with
GP582R & GP586R

GP587SU



M Transversal
to be used with
GP582R & GP586R

GP588SU



M Transversal
to be used with
GP582R & GP586R

GP591SU



L Transversal
to be used with
GP590R

GP594SU




L Transversal
to be used with
GP590R

GP595SU



L Transversal
to be used with
GP590R

GP598SU



L Transversal
to be used with
GP590R



To be used with

GP580R


65 mm



S Transversal Saw Blades

GP582R


65 mm



M Transversal Saw Blades

GP586R

95 mm



M Transversal Saw Blades





GP590R

30 mm




L Transversal Saw Blades

Burrs



ELAN[®] 4 Power System

Burrs - 1 Ring



Ø	0.6 mm	0.8 mm	1.0 mm	1.0 mm	1.4 mm	1.4 mm	1.8 mm	1.8 mm	2.3 mm	2.3 mm	2.7 mm	
I	GP111SU	GP112SU	GP113SU	GP114SU	GP115SU	GP116SU	GP117SU	GP118SU	GP120SU	GP121SU	GP122SU	Rosen


Ø	2.7 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm	
I	GP123SU	GP124SU	GP125SU	GP126SU	GP127SU	GP128SU	GP129SU	Rosen

Ø	3.0 mm	4.0 mm	5.0 mm	6.0 mm	
Soft Cut 					
I	GP133SU	GP134SU	GP135SU	GP136SU	Rosen

Ø	4.0 mm	5.0 mm	6.0 mm	7.0 mm	
I	GP184SU	GP185SU	GP186SU	GP187SU	Twin-Cut Burr



Ø	0.6 mm	0.8 mm	1.0 mm	1.0 mm	1.4 mm	1.4 mm	1.8 mm	1.8 mm	2.3 mm	2.3 mm	2.7 mm	
I	GP141SU	GP142SU	GP143SU	GP144SU	GP145SU	GP146SU	GP147SU	GP148SU	GP149SU	GP150SU	GP151SU	Diamond Burr

Ø	2.7 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm	
I	GP152SU	GP153SU	GP154SU	GP155SU	GP156SU	GP157SU	GP158SU	Diamond Burr



Ø	2.0 mm	2.0 mm	2.3 mm	2.3 mm	3.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm	
										 sterile
I	GP161SU	GP162SU	GP163SU	GP164SU	GP165SU	GP166SU	GP168SU	GP169SU	GP170SU	Diamond Burr, Coarse

Ø	3.0 mm	4.0 mm	5.0 mm	6.0 mm	7.0 mm	
						 sterile
I	GP173SU	GP174SU	GP175SU	GP176SU	GP177SU	Diamond Burr, Extra Coarse

Ø	1.5 mm	2.0 mm	2.5 mm	3.0 mm	3.0 mm
US Type 					
I	GP201SU	GP202SU	GP203SU	GP204SU	GP205SU

 sterile
Neuro Cutter

Ø	2.0 mm	2.5 mm	3.0 mm
Aesculap Type 			
I	GP208SU	GP209SU	GP210SU


 sterile
Neuro Cutter

Ø	1.5 mm	1.5 mm	2.0 mm	2.0 mm	2.5 mm	2.5 mm	3.0 mm	3.0 mm	4.0 mm	4.0 mm	
											 sterile
I	GP211SU	GP212SU	GP213SU	GP214SU	GP215SU	GP216SU	GP217SU	GP218SU	GP219SU	GP220SU	Neuro Cutter Diamond



ELAN[®] 4 Power System

Burrs - 1 Ring



Ø	4.0 mm	5.0 mm	6.0 mm
Standard			
	GP224SU	GP225SU	GP226SU

sterile
Barrel Burr

Ø	4.0 mm	5.0 mm	6.0 mm
Soft Cut			
	GP228SU	GP229SU	GP230SU

sterile
Barrel Burr

Ø	4.0 mm	5.0 mm	6.0 mm
	GP234SU	GP235SU	GP236SU

sterile
Cone Burr

Ø	3.0 mm	4.0 mm	5.0 mm	6.0 mm
Midas Type				
	GP193SU	GP194SU	GP195SU	GP196SU

sterile
Acorn Burr

Ø	4.0 mm	5.5 mm
	GP232SU	GP233SU

sterile
Oval Burr

Ø	4.0 mm
	GP238SU

sterile
Reverse Taper Burr Coarse Diamond

Ø	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm
Tungsten Carbide						
	GP113TC	GP119TC	GP124TC	GP126TC	GP128TC	GP129TC

sterile
Rosen

Ø	3.0 mm
Tungsten Carbide	
	GP210TC

sterile
Neuro Cutter

Ø	2.3 mm	1.8 mm	2.1 mm	2.3 mm
	GP240SU	GP241SU	GP242SU	GP243SU


sterile
Lindemann

Ø	1.5 mm	2.0 mm
	GP188SU	GP189SU

sterile
Twist Drill


Ø	1.0 mm
	GP190SU

sterile
Pin Cutter




Burrs – 2 Ring

Ø	2.3 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm	7.0 mm	8.0 mm	
										 sterile
II	GP301SU	GP302SU	GP303SU	GP304SU	GP305SU	GP306SU	GP307SU	GP308SU	GP309SU	Rosen

Ø	2.3 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm	
								 sterile
II	GP311SU	GP312SU	GP313SU	GP314SU	GP315SU	GP316SU	GP317SU	Diamond Burr

Ø	2.3 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm	
						 sterile
II	GP321SU	GP322SU	GP323SU	GP324SU	GP325SU	Diamond Burr, Coarse

Ø	4.0 mm	5.0 mm	6.0 mm	
				 sterile
II	GP328SU	GP329SU	GP330SU	Diamond Burr, Extra Coarse

Ø	2.0 mm	2.5 mm	3.0 mm	3.0 mm
US Type				
II	GP361SU	GP362SU	GP363SU	GP364SU


 sterile

Neuro Cutter

Ø	3.0 mm
Aesculap Type	
II	GP367SU

 sterile

Neuro Cutter

Ø	3.0 mm	3.0 mm
		
II	GP370SU	GP371SU



 sterile

Neuro Cutter Diamond




ELAN[®] 4 Power System


Burrs – 2 Ring




Barrel Burr




Barrel Burr




Cone Burr




Acorn Burr




Oval Burr




Reverse Taper Burr Coarse Diamond




Craniotome Cutter – Spiral Type




Craniotome Cutter – Straight Type




Lindemann



Twist Drill



Twist Drill



Ø	1.0 mm	Ø	5.0 mm	6.0 mm	7.0 mm	8.0 mm	9.0 mm
							
II	GP350SU	II	GP335SU	GP336SU	GP337SU	GP338SU	GP339SU
 sterile		 sterile					
Pin Cutter		Twin-Cut Burr					

Tungsten Carbide Burrs

Ø	2.3 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm	7.0 mm	8.0 mm
							
II	GP301TC	GP302TC	GP304TC	GP306TC	GP307TC	GP308TC	GP309TC
 sterile							
Rosen							

Ø	1.2 mm	1.6 mm	Ø	3.0 mm
				
II	GP395TC	GP396TC	II	GP367TC
 sterile			 sterile	
Taper Burr			Neuro Cutter	




Burrs for metal dissection

Ø	2.3 mm
	
II ②	GP399TC-SU
 sterile	
Side Cutter	





ELAN[®] 4 Power System

MIS Burrs












10 cm















13 cm



Ø	3.0 mm	4.0 mm	Ø	3.0 mm	4.0 mm	Ø	2.0 mm	3.0 mm	Ø	2.0 mm	3.0 mm
											
L10	GP422SU	GP423SU	L13	GP442SU	GP443SU	L10	GP420SU	GP421SU	L13	GP440SU	GP441SU
 sterile			 sterile			 sterile			 sterile		
Rosen Burr			Rosen Burr			Neuro Cutter			Neuro Cutter		

Ø	2.0 mm	3.0 mm	4.0 mm	Ø	2.0 mm	3.0 mm	4.0 mm	5.0 mm
								
L10	GP424SU	GP425SU	GP426SU	L13	GP444SU	GP445SU	GP446SU	GP447SU
 sterile				 sterile				
Neuro Cutter Diamond				Neuro Cutter Diamond				

Ø	3.0 mm	4.0 mm	5.0 mm	Ø	3.0 mm	4.0 mm	5.0 mm	Ø	3.0 mm	4.0 mm	5.0 mm
											
L10	GP432SU	GP433SU	GP434SU	L13	GP452SU	GP453SU	GP454SU	L10	GP435SU	GP436SU	GP437SU
 sterile				 sterile				 sterile			
Diamond Burr				Diamond Burr				Diamond Burr Extra Coarse			

Ø	3.0 mm	4.0 mm	5.0 mm	Ø	3.0 mm	4.0 mm	5.0 mm	Ø	3.0 mm	4.0 mm	5.0 mm
											
L13	GP455SU	GP456SU	GP457SU	L10	GP428SU	GP429SU	GP430SU	L13	GP448SU	GP449SU	GP450SU
 sterile				 sterile				 sterile			
Diamond Burr Extra Coarse				Neuro Cutter Diamond Extra Coarse				Neuro Cutter Diamond Extra Coarse			





Through collaborative
excellence we will improve the
quality of a patient's life and
meet the needs of the changing
healthcare environment.

Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form.
In the case of misuse we retain the rights to recall our catalogs and price lists and to take legal actions.

©2023 AESCULAP, INC. ALL RIGHTS RESERVED. PRINTED IN THE USA.
Aesculap is an equal opportunity employer

Aesculap, Inc. | 3773 Corporate Parkway | Center Valley, PA | 18034
Phone 800-282-9000 | Fax 610-791-6886 | www.aesculapusa.com

Aesculap, Inc. – a B. Braun company

DOC1360 Rev F 082323

