



Sterile Container Systems

General Catalog

OPERATE WITH GREATER PRECISION™



AESCULAP®

Table of Contents

Introduction	4-5
AESCU LAP Aicon® Sterile Container System Containers and Baskets.....	6-17
▪ AEscuLAP Aicon Sterile Container System Design	6-9
▪ Full-Size	10-11
▪ Three-Quarter Size	12-13
▪ Half-Size	14-15
▪ Identification Faceplates and Tags	16
▪ Processing Supplies and Accessories	17
SterilContainer® System Containers and Baskets	18-37
▪ SterilContainer Design, , IUSS Validated	18-19
▪ Full-Size	20-21
▪ Three-Quarter Size	22-23
▪ Half-Size	24-25
▪ Extra-Long Mini Size	26
▪ Extra-Long Size	27
▪ Quarter Size, , Bar Code ID Brackets for Baskets	28-29
▪ Mini Size	30
▪ Wide-Body Size	31
▪ JS Series S2 and JM SterilContainer™ S Container Systems	32-33
▪ Sterile Processing Supplies	34-36
▪ Identification Labels and Tags	37
Instrument Organization System	38-51
General Accessories	52-61
▪ Endoscopy Storage Accessories	52-53
▪ Instrument Storage Racks	54-55
▪ Basket Lid Instrument Support Holder	56
▪ Micro Instruments Holders and Racks	57
▪ Mats, Lubricants and Cleaners	58-59
▪ Communicator Tags, Stringers and ID Brackets	60
Storage and Transport	61
▪ Display Peg Boards and Instrument Holders	61
Miscellaneous	62-64
▪ Metal Basins and Bowls	62
▪ SterilContainer Replacement Parts	63
▪ BIO-BIN Containers	64
Sterile Container Product Dimensions	66-79
▪ AEscuLAP Aicon Sterile Container System	66-69
▪ SterilContainer System	70-79
Instructions For Use	80

Commitment to Quality



A Heritage of Quality for 150 Years



Gottfried Jetter

Aesculap, the world's largest manufacturer of surgical instruments, has been crafting high performance instruments in Tuttlingen, Germany since 1867. Named

for the Greek god of healing, Asklepios, the company was founded by Gottfried Jetter, the first person ever to produce standardized surgical instruments.

Today, Aesculap employs nearly 12,000 people in our modern, Aesculap-owned and operated factories around the world, including 3,500 people in the Tuttlingen facilities which house over 450,000 square feet of production area.

PERFORMANCE QUALITY Surgical Instruments Made by Aesculap

For every surgical instrument *Made by Aesculap* we demand not just "Quality", but *PERFORMANCE QUALITY*. A surgical instrument with PERFORMANCE QUALITY is meticulously crafted to have a solid feel and to perform like an extension of the surgeon's own hand. To deliver PERFORMANCE QUALITY worthy of the Aesculap name, we start with the finest, raw materials including highly pure stainless steel which is purchased and tested in Tuttlingen. Our master craftsmen, who are among the most rigorously trained and experienced in the industry, then follow the strictest standards when crafting each instrument, including:

- **State-of-the-art German-made equipment** which is used in every Aesculap-owned factory worldwide
- **Aesculap hand-crafted, Gold-Plated Master QC Standards** that assure pattern consistency
- **Strict Aesculap manufacturing guidelines**, including tight tolerances more strict than International Quality Standards, for precision and reliability from first use and throughout the lifetime of each instrument

As our Mission states, we are "passionately committed to providing high quality, innovative products and services to all surgical disciplines."

Aesculap Product Range:

- Surgical instruments for all disciplines
- Sterile containers
- Endoscopy
- Neuroendoscopy
- Aneurysm clips
- Power systems
- Pneumatic KAIRisons
- Spinal implants
- Orthopaedics
- Per-procedure programs
- Repair services
- Consulting services
- Dura substitutes
- Hydrocephalus

Aesculap U.S. Headquarters



Aesculap U.S. Headquarters

Our U.S. Headquarters, located in Center Valley, Pennsylvania, was established in 1977 to serve the needs of the U.S. market.

Should you have any need or questions, our customer service professionals welcome your call.

Aesculap holds customers and patients to the highest level in our core values. Our commitment is to exceed your expectations.

Please Call Customer Service For:

- Placing orders
- Product information
- Pricing

Aesculap, Inc.
3773 Corporate Parkway
Center Valley, PA 18034
Phone 800-258-1946
Fax 610-791-6880
www.aesculapusa.com

Customer Service:

Phone: 800-282-9000

Fax: 610-791-6886

AESCLAP Aicon®

Sterile Container System

USER-FRIENDLY CLOSURE SYSTEM

Single-motion latch with large lifting handle designed to improve lid handling during aseptic presentation.

STERILE SET SECURITY

Tamper-evident lock is visible from the front and top of the container and lays flat to protect it during transportation and handling.

EASY ACCESS TO INFORMATION


QR code provides product information and Instructions for Use (IFU) by simply scanning the code with a mobile device camera.

IMPROVED SET IDENTIFICATION

Interchangeable color faceplates and tags allow for more than 100 set identification options designed to easily identify specialties. Retaining brackets secure additional information.

ERGONOMIC CONTAINER HANDLING

Integrated handle evenly distributes weight through whole hand to help reduce fatigue and designed for greater access to rear handle during lifting.





reddot winner 2020

ADDED STERILITY PROTECTION

Optional cover designed for additional protection when transporting between facilities and may remain in place during cleaning and sterilization.

REDUCED PROCESSING SUPPLIES


Full-sized container only requires one filter per set, helping to reduce processing supply expenses. Filter retention plate holds filter in place to help protect filter from damage and strike through.

LONGER SERVICE LIFE

Robust container and gasket designed to extend container life and provide consistent performance.

ENHANCED DRYING SYSTEM (EDS)

Models with EDS are designed to reduce set drying times, decrease overall set reprocessing time and provide additional confidence of delivering dry sets to the operating room.




AESCULAP Aicon®

Sterile Container System

INSTRUMENT SECURITY

Non-slipping wavelike basket bottom designed to hold instruments in place. Ring instruments stay upright during the packaging and transportation process to help protect instruments and make them more accessible.



ADDED PROTECTION

Optional basket lid designed for additional protection of set contents.

ERGONOMIC BASKET HANDLING


Large handles designed to improve user control in the Sterile Processing Department and during aseptic presentation in the Operating Room.


QUICK SET IDENTIFICATION

Numerous tag options designed to allow for easy set labeling and quick set identification.

BASKET FLEXIBILITY


Smooth and rounded design allows basket to be used with containers or single-use wrap.






BASKET AND CONTAINER SYNCHRONIZATION

Enables identical tags to be used on the basket and removable container faceplate to streamline the reprocessing cycle and help improve set identification.



SPACE-SAVING BASKET DESIGN

Baskets nest together securely, with or without basket lid, to allow stacking in containers and during storage.



AESCULAP Aicon®

Full-Size Sterile Container System

CONTAINER BOTTOMS

Item no.	Length Width	23 ⁵ / ₈ in 11 ³ / ₄ in	598 mm 297 mm	Height ¹
JJ110				4 in 100 mm
JJ120				6 in 150 mm
JJ130				8 in 200 mm
JJ140				10 in 250 mm

CONTAINER BOTTOMS WITH ENHANCED DRYING SYSTEM (EDS)



Item no.	Length Width	23 ⁵ / ₈ in 11 ³ / ₄ in	598 mm 297 mm	Height ¹
JJ111				4 in 100 mm
JJ121				6 in 150 mm
JJ131				8 in 200 mm
JJ141				10 in 250 mm

LIDS FIT ALL FULL-SIZE CONTAINER BOTTOMS

Item no.	Length Width	23 ⁵ / ₈ in 11 ³ / ₄ in	598 mm 297 mm	Height
JJ410				1 ¹ / ₈ in 28 mm

¹ Height with lid. Inches rounded up to closest 1/8 inch.

AESCULAP Aicon®

Full-Size Sterile Container System Basket Sizes



BASKET COMBINATION OPTIONS

Item no. container bottom	Item no. baskets FULL-SIZE
JJ110 JJ111	
JJ120 JJ121	
JJ130 JJ131	
JJ140 JJ141	

Container bottom
 Basket
 Front view

FULL-SIZE BASKET SIZES

Item no.	Length 21¼ in 537 mm Width 10 in 250 mm	Height¹
JJ161R	1/1 AESCULAP Aicon basket	1½ in 35 mm
JJ162R		2¼ in 55 mm
JJ163R		3¼ in 80 mm
JJ164R		4¼ in 105 mm
JJ167R	1/1 AESCULAP Aicon basket lid	¾ in 18 mm

¹ Inches rounded up to closest 1/8 inch.

AESCULAP Aicon®

Three-Quarter Size Sterile Container System

CONTAINER BOTTOMS

Item no.	Length 18 ³ / ₈ in 464 mm Width 11 ³ / ₄ in 297 mm	Height ¹
JJ310		4 in 100 mm
JJ320		6 in 150 mm
JJ330		8 in 200 mm
JJ340		10 in 250 mm

CONTAINER BOTTOMS WITH ENHANCED DRYING SYSTEM (EDS)



Item no.	Length 18 ³ / ₈ in 464 mm Width 11 ³ / ₄ in 297 mm	Height ¹
JJ311		4 in 100 mm
JJ321		6 in 150 mm
JJ331		8 in 200 mm
JJ341		10 in 250 mm

LIDS FIT ALL THREE-QUARTER SIZE CONTAINER BOTTOMS

Item no.	Length 18 ³ / ₈ in 464 mm Width 11 ³ / ₄ in 297 mm	Height
JJ430		1 ¹ / ₈ in 28 mm

¹ Height with lid. Inches rounded up to closest 1/8 inch.

AESCULAP Aicon®

Three-Quarter Size Sterile Container System Basket Sizes

BASKET COMBINATION OPTIONS



Item no. container bottoms	Item no. baskets THREE-QUARTER SIZE
JJ310 JJ311	
JJ320 JJ321	
JJ330 JJ331	
JJ340 JJ341	

Container bottom
 Basket
 Front view

THREE-QUARTER SIZE BASKET SIZES

Item no.	Length 15 ⁷ / ₈ in 403 mm Width 10 in 250 mm	Height ¹
JJ361R	3/4 AESCULAP Aicon basket	1 ¹ / ₂ in 35 mm
JJ362R		2 ¹ / ₄ in 55 mm
JJ363R		3 ¹ / ₄ in 80 mm
JJ364R		4 ¹ / ₄ in 105 mm
JJ367R	3/4 AESCULAP Aicon basket lid	3/4 in 18 mm

¹ Inches rounded up to closest 1/8 inch.

AESCULAP Aicon®

Half-Size Sterile Container System

CONTAINER BOTTOMS

Item no.	Length 11 ⁷ / ₈ in 300 mm Width 11 ³ / ₄ in 297 mm	Height ¹
JJ210		4 in 100 mm
JJ220		6 in 150 mm
JJ230		8 in 200 mm
JJ240		10 in 250 mm

CONTAINER BOTTOMS WITH ENHANCED DRYING SYSTEM (EDS)



Item no.	Length 11 ⁵ / ₈ in 297 mm Width 11 ³ / ₄ in 300 mm	Height ¹
JJ211		4 in 100 mm
JJ221		6 in 150 mm
JJ231		8 in 200 mm
JJ241		10 in 250 mm

LIDS FIT ALL HALF-SIZE CONTAINER BOTTOMS



Item no.	Length 11 ⁷ / ₈ in 300 mm Width 11 ³ / ₄ in 297 mm	Height
JJ420		1 ¹ / ₈ in 28 mm

¹ Height with lid. Inches rounded up to closest 1/8 inch.

AESCULAP Aicon®

Half-Size Sterile Container System Basket Sizes

BASKET COMBINATION OPTIONS

Item no. container bottoms	Item no. baskets HALF-SIZE	
JJ210 JJ211	JJ262R	
JJ220 JJ221	JJ264R, JJ261R, JJ262R, JJ261R	
JJ230 JJ231	JJ262R, JJ263R, JJ264R, JJ263R	
JJ240 JJ241	JJ264R	

 Container bottom
  Basket
 Front view

HALF-SIZE BASKET SIZES

Item no.	Length Width	9 1/2 in 10 in	240 mm 250 mm	Height ¹
JJ261R	1/2 AESCULAP Aicon basket			1 1/2 in 35 mm
JJ262R				2 1/4 in 55 mm
JJ263R				3 1/4 in 80 mm
JJ264R				4 1/4 in 105 mm
JJ267R	1/2 AESCULAP Aicon basket lid			3/4 in 18 mm

¹ Inches rounded up to closest 1/8 inch.

AESCLAP Aicon®

Sterile Container System Identification Faceplates and Tags

Faceplates required if SPD utilizing Identification Cards.

JJ582 Adaptor designed for use with Aesculap Indicator Cards MD346/MD345 located on page 34.

Color		Faceplate		Tag	
		1/pk	10/pk	Large, 1/pk	Small, 1/pk
Red		JJ500	JJ520	JJ540U	JJ560U
Blue		JJ501	JJ521	JJ541U	JJ561U
Green		JJ502	JJ522	JJ542U	JJ562U
Yellow		JJ503	JJ523	JJ543U	JJ563U
Gray		JJ504	JJ524	JJ544U	JJ564U
Black		JJ505	JJ525	-	-
White		JJ506	JJ526	JJ546U	JJ566U
Orange		JJ507	JJ527	JJ547U	JJ567U
Violet		JJ508	JJ528	JJ548U	JJ568U
Brown		JJ509	JJ529	JJ549U	JJ569U
Size		9 7/8 x 2 3/4 in 248 x 70 mm		2 1/4 x 1 1/4 in 56 x 29 mm	2 1/4 x 3/4 in 56 x 19 mm



Faceplates



Large Tags

2 1/4 x 1 1/4 in
56 x 29 mm



Small Tags


2 1/4 x 3/4 in
56 x 19 mm

Bar Code ID Brackets for Baskets



JJ586U

96 x 34 mm (3 3/4" x 1 1/4")



JJ587U

105 x 54 mm (4 1/8" x 2 1/8")



AESCLAP Aicon®

Sterile Container System Processing Supplies and Accessories

SINGLE-USE PAPER FILTER

Item no. 100/pk	Without indicator
JJ612	

RETENTION PLATE

Item no. 1/pk	
JJ600	

COVER FOR PERFORATED AREA

Item no. 1/pk	Item no. 10/pk	
JJ440	JJ441	

Cover designed for transport of used and soiled instruments.

*Refer to AESCLAP Aicon IFU for further information. www.aesculapusaifus.com

RETAINING BRACKET FOR FACEPLATE

Item no. 10/pk	
JJ599	

LABEL ADAPTER FOR FACEPLATE

Item no. 1/pk	3 1/4 x 1 1/2 in 80 x 35 mm
JJ582	

Used in conjunction with Aesculap Indicator Cards MD345/MD346 Pg. 34

CONTAINER LOCK

Item no. 100/pk	Item no. 1000/pk	Blue, without indicator
JJ705	JJ700	
Item no. 100/pk	Item no. 1000/pk	Orange, with steam indicator
JJ725	JJ720	
Item no. 100/pk	Item no. 1000/pk	Yellow, with steam indicator
JJ726	JJ721	

ENHANCED DRYING SYSTEM (EDS) MODULE

Item no. 1/pk	
JJ800	

BIOHAZARD CARD

Item no. 250/box	
US991	

For use with soiled instruments.

*EDS Aicon version not approved for BioBin.

SterilContainer™

System



Validated for
IUSS

Contact your
Sales Representative
for information

BENEFITS

- Extensive range of container sizes and heights designed to meet all your instrument sets requirements
- Lightweight aluminum construction helps you meet AAMI 25lbs weight recommendation
- 90° Stop handle designed to protect staff from finger injuries during handling
- Single motion latch designed to secure containers quickly to help reduce time and hassle in the SPD and OR
- One step retention plate with audible "click" designed to secure filter to simplify operations and help improve processing time
- Able to containerize pre-configured instrument sets including orthopaedic, spine and specialty medical devices
- Solid bottom configured with either Lid Option "VALIDATED for IUSS"

As pictured - Validation with Primeline Pro Lid Reusable Filter for IUSS



SterilContainer™

Full-Size System

Full-Size Containers – Length 23¼ x width 11¼

Perforated Bottoms with Retention Plates



JN440
4 ¼-inch height



JN441
5 ½-inch height



JN442
6-inch height



JN444
8-inch height



JN446
10 ½-inch height

*Basket height includes feet

Solid Bottoms



JK440



JK441



JK442




JK444



JK446

Full-Size Container Lids

Aluminum SterilContainer™ Lid



- JK485 Full-size red lid
- JK486 Full-size blue lid
- JK487 Full-size green lid
- JK488 Full-size gold lid
- JK489 Full-size silver lid

Note: Filter requirements when used with aluminum lids:
Perforated full-size container with lids require *four* round filters**.
Solid full size container with lids require *two* round filters**.

PrimeLine® Pro Full-Size Container Lids



- JP101 Full-size red lid
- JP102 Full-size blue lid
- JP103 Full-size green lid
- JP104 Full-size gold lid
- JP105 Full-size silver lid

Note: Integrated Reusable PTFE Filter (JP050 included with Lid Purchase) designed for 2200 sterilization cycles.

*All dimensions shown here are outside dimensions with the lid.

SterilContainer™

Full-Size System Basket Sizes and Accessories

Full-Size Stainless Steel Basket Standard Perforation



Length x Width x Height (in.)
Outside dimensions

JF221R 21¼ x 10 x 1½*

JF222R 21¼ x 10 x 2¼*

JF522R 21¼ x 10 x 2¼**

JF223R 21¼ x 10 x 3*

JF224R 21¼ x 10 x 4¼*

*Basket height includes feet

**Telescoping with 2¾" handle

Basket Cover for Perforated Baskets



Length x Width x Height (in.)
Outside dimensions

JF227R 2½ x 10 1/16 x ¾

Perforated Basket Cover
for JF Baskets


FULL-SIZE BASKET SIZES

Item no.	Length Width	21¼ in 10 in	537 mm 250 mm	Height ¹
JJ161R	1/1 AESCULAP Aicon® basket			1½ in 35 mm
JJ162R				2¼ in 55 mm
JJ163R				3¼ in 80 mm
JJ164R				4¼ in 105 mm
JJ167R	1/1 AESCULAP Aicon basket lid			¾ in 18 mm



Silicone Mesh Mats & Silicone Holding Mats with Perforations

Silicone Mesh Mats





JF938 - 9½" x 9⅞" Half-size silicone mat

JF939 - 18⅞" x 9⅞" Din silicone mat

JF940 - 15⅞" x 9⅞" Three-quarter size silicone mat

JF941 - 21⅞" x 9⅞" Full-size silicone mat

Perforated Silicone Holding Mats



JF944 - 9¾" x 4" Mini-size silicone mat

JF945 - 9¾" x 9¼" Half-size silicone mat

JF947 - 15¼" x 9¼" Three-quarter size silicone mat

JF949 - 20¼" x 9¼" Full-size silicone mat

*All dimensions shown here are outside dimensions with the lid.

¹ Inches rounded up to closest 1/8 inch. Detailed dimensions on page 14-15.

SterilContainer™

Three-Quarter Size System

Three-Quarter Size Containers – Length 18½ x width 11¼

Perforated Bottoms with Retention Plates



JN740
4¼-inch height



JN741
5½-inch height



JN742
6-inch height



JN744
8-inch height

Solid Bottoms



Three-Quarter Size Container Lids

Aluminum SterilContainer™ Lid



- JK785 Three-quarter size red lid
- JK786 Three-quarter size blue lid
- JK787 Three-quarter size green lid
- JK788 Three-quarter size gold lid
- JK789 Three-quarter size silver lid

Note: Filter requirements when used with aluminum lids:
Perforated three-quarter size container with lids require *two* round filters**.
Solid three-quarter size container with lids require *one* round filter**.

PrimeLine® Pro Three-Quarter Size Container Lids



- JP111 Three-quarter size red lid
- JP112 Three-quarter size blue lid
- JP113 Three-quarter size green lid
- JP114 Three-quarter size gold lid
- JP115 Three-quarter size silver lid

Note: Integrated Reusable PTFE Filter (JP050 included with Lid Purchase) designed for 2200 sterilization cycles.

*All dimensions shown here are outside dimensions with the lid.

SterilContainer™

Three-Quarter System Basket Sizes and Accessories

Three-Quarter Size Stainless Steel Basket Standard Perforation



Length x Width x Height (in.)
Outside dimensions

JF251R 16 x 10 x 1½*

JF252R 16 x 10 x 2¼*

JF253R 16 x 10 x 3*

JF254R 16 x 10 x 4⅛*

*Basket height includes feet

Perforated Basket

- Fold-in handles
- Plastic feet

Basket Cover for Perforated Baskets



Length x Width x Height (in.)
Outside dimensions

JF257R 16½ x 10⅞ x ¾


Perforated Basket Cover for JF Baskets

THREE-QUARTER BASKET SIZES

Item no.	Length Width	15⅞ in 10 in	403 mm 250 mm	Height¹
JJ361R	3/4 AESCULAP Aicon® basket			1½ in 35 mm
JJ362R				2¼ in 55 mm
JJ363R				3¼ in 80 mm
JJ364R				4¼ in 105 mm
JJ367R	3/4 AESCULAP Aicon basket lid			¾ in 18 mm



Silicone Mesh Mats & Silicone Holding Mats with Perforations





JF938 – 9½" x 9⅞" Half-size silicone mat

JF939 – 18⅞" x 9⅞" Din silicone mat

JF940 – 15⅞" x 9⅞" Three-quarter size silicone mat

JF941 – 21⅞" x 9⅞" Full-size silicone mat

Perforated Silicone Holding Mats



JF944 – 9¾" x 4" Mini-size silicone mat

JF945 – 9¾" x 9¼" Half-size silicone mat

JF947 – 15¼" x 9¼" Three-quarter size silicone mat

JF949 – 20¼" x 9¼" Full-size silicone mat

*All dimensions shown here are outside dimensions with the lid.

¹Inches rounded up to closest 1/8 inch. Detailed dimensions on page 14–15.

SterilContainer™

Half-Size System

Half-Size Containers – Length 11¾ x width 11¼

Perforated Bottoms with Retention Plates



JN340

4¼-inch height



JN341

5½-inch height



JN342

6-inch height



JN344

8-inch height



JN346

10½-inch height

Solid Bottoms



JK340



JK341



JK342




JK344



JK346

Half-Size Container Lids

Aluminum SterilContainer™ 2000 Lid



- JK385 Half-size red lid
- JK386 Half-size blue lid
- JK387 Half-size green lid
- JK388 Half-size gold lid
- JK389 Half-size silver lid

Note: Filter requirements when used with aluminum lids:
Perforated half-size container with lids require *two* round filters**.
Solid half-size container with lids require *one* round filter**.

PrimeLine® Pro Half- Size Container Lids



- JP121 Half-size red lid
- JP122 Half-size blue lid
- JP123 Half-size green lid
- JP124 Half-size gold lid
- JP125 Half-size silver lid



Note: Integrated Reusable PTFE Filter (JP050 included with Lid Purchase) designed for 2200 sterilization cycles.

*All dimensions shown here are outside dimensions with the lid.

SterilContainer™

Half-Size System Basket Sizes and Accessories

Half-Size Stainless Steel Basket Standard Perforation



JF112R

Perforated Basket

- Fold-in handles
- Plastic feet

Length x Width x Height (in.)
Outside dimensions

JF111R 10 x 9⁷/₈ x 1¹/₂*


JF112R 10 x 9⁷/₈ x 2¹/₄*

JF113R 10 x 9⁷/₈ x 3*

JF114R 10 x 9⁷/₈ x 4¹/₈*

*Basket height includes feet

Basket Cover for Perforated Baskets



JF117R

Perforated Basket Cover for JF Baskets

Length x Width x Height (in.)
Outside dimensions


JF117R 10¹/₈ x 9³/₄ x ³/₄

HALF-SIZE BASKET SIZES

Item no.	Length Width	9 ¹ / ₂ in 10 in	240 mm 250 mm	Height ¹
JJ261R	1/2 AESCULAP Aicon® basket			1 ¹ / ₂ in 35 mm
JJ262R				2 ¹ / ₄ in 55 mm
JJ263R				3 ¹ / ₄ in 80 mm
JJ264R				4 ¹ / ₄ in 105 mm
JJ267R	1/2 AESCULAP Aicon basket lid			³ / ₄ in 18 mm



Silicone Mesh Mats & Silicone Holding Mats with Perforations





JF938 – 9¹/₂" x 9⁷/₈" Half-size silicone mat

JF939 – 18⁷/₈" x 9⁷/₈" Din silicone mat

JF940 – 15⁷/₈" x 9⁷/₈" Three-quarter size silicone mat

JF941 – 21¹/₈" x 9⁷/₈" Full-size silicone mat

Perforated Silicone Holding Mats



JF944 – 9³/₄" x 4" Mini-size silicone mat

JF945 – 9³/₄" x 9¹/₄" Half-size silicone mat

JF947 – 15¹/₄" x 9¹/₄" Three-quarter size silicone mat

JF949 – 20¹/₄" x 9¹/₄" Full-size silicone mat

*All dimensions shown here are outside dimensions with the lid.

¹ Inches rounded up to closest 1/8 inch. Detailed dimensions on page 14–15.

SterilContainer™


Extra-Long Mini Size Container System

Extra-Long Mini Size Containers – Length 20 x Width 5½

Non-Anodized



JM021
3-inch height



JM020 Lid




MD150 Lid with Handle

Anodized for Steam and EtO



JN021
3-inch height




JK020 Lid




MD153 Lid with Handle

Extra-Long Mini Size Container Baskets



JF435R
Basket Only



JF428R
Single Scope



JF429R
Dual Scope

Length x Width x Height (in.) Outside dimensions

JF435R 20 x 5¼ x 1¾

JF428R 20 x 5¼ x 1¾

JF429R 20 x 5¼ x 1¾

Extra-Long Mini Size Container Accessories



JG367
Silicone Scope Eyepiece Holder



JG368
Silicone Scope Sleeve Holder



JG369
Support for Angled Optics



JG340
Small Silicone Holder for 5 mm instruments



JG343
Small Silicone Holder for 10 mm instruments

JG300/JG303R fixation screws sold separately

*All dimensions shown here are outside dimensions with the lid.

SterilContainer™

Extra-Long Container System

Extra-Long Containers – Length 26¼ x Width 11¼

Perforated Bottom



JN445

8-inch height

Solid Bottom



JK443

5½-inch height

Perforated Bottom



JN443

5½-inch height

Lid with Retention Plates



JK490

JK490 – Extra-Long Container Lid

Extra-Long Perforated Basket

Perforated Basket



JF232R

Length x Width x Height (in.)
Outside dimensions


JF232R 26 x 10 x 2¼*

JF532R 26 x 10 x 2¼ with 2¾" telescoping handle*

JF236R 26 x 10 x 3⅜ will not work with PL964R rack*


*Basket height includes feet

Extra-Long Container Accessories




JF933

Extra-Long Silicone Mat



MD375

Endorack – Extra Long,
Holds 14 Instruments.
Fits in JN445 8" Container




MD379

Endorack – Extra Long Kit
Includes MD375 Rack, JF933 Mat,
JF232R Basket



PL961R

PL961R
Endorack Extra-Long,
perforated basket and
silicone mesh mat, for
instruments with shaft-length up to 42
cm, fits into JN445



PL964R

PL964R
Endorack Extra-Long, fits into 8" height,
Extra-Long SterilContainer

*All dimensions shown here are outside dimensions with the lid.

SterilContainer™

Quarter-Size System

Quarter-Size Containers – Length 12½ x Width 7½ (Bottom and Lid Sold Separately)

Container System

(Perforated Bottoms Only)

Baskets



JN086

1⅞-inch height



JG381R

Quarter-Size Perforated Basket with Lid and Three Instrument Racks
Length x Width x Height – Outside (in.)
10¾ x 7 x 1⅜



JN088

2½-inch height



JF264R

Quarter-Size Perforated Basket with Lid
Length x Width x Height – Outside (in.)
10¾ x 7 x 1⅝



MD374

Quarter-Size Perforated Basket with Lid and Silicone Cushioning Rack
Length x Width x Height – Outside (in.)
10¾ x 7 x 1⅝



JG384R

Quarter-Size Perforated Basket with Lid and Three Instrument Racks
One JG381R and one JG384R can fit into JN095
Length x Width x Height – Outside (in.)
10¾ x 7 x 1⅝



JN089

3⅝-inch height



JG285R

Quarter-Size Perforated Basket with Lid
One can fit into JN095
Length x Width x Height – Outside (in.)
10¾ x 6¾ x 2⅜



JN090

5⅛-inch height



JF282R

Quarter-Size Tray
fits into container JN096
Length x Width x Height – Outside (in.)
10⅝ x 6¾ x 3⅝



JF289R

Cover for JF282R

Quarter-Size Lid for tray JF282R
Length x Width x Height – Outside (in.)
10¾ x 7 x ⅝

*All dimensions shown here are outside dimensions with the lid.

Filter Requirements: US999 or JK091 filters.

SterilContainer™

Quarter-Size System Lid Sizes and Accessories

Quarter-Size Aluminum Container Lids




with Handles

Note: MD151 does not include frames for labels.

Bar Code ID Brackets for Baskets

JG616U and JG617U are Bar Code Identification brackets that can be attached to AESCULAP Aicon®, JJ Series, JF Series or perforated baskets to aid in securing bar code labels to the basket. The stainless steel brackets are secured to the baskets by bending the flexible tabs into the perforations.



*All dimensions shown here are outside dimensions with the lid.

SterilContainer™

Mini-Size System

Mini-Size Containers – Length 12¼ x Width 5½

Perforated Bottoms with Retention Plate



JN187

1½-inch height

NOTE: Shallow Mini container does not have an internal basket, however a Mini mat can be used.



JN188

2⅝-inch height

Solid Bottoms



JK187

1½-inch height




JK188

2⅝-inch height

Mini-Size Container Lids

Aluminum SterilContainer™ Lid



JK170 Mini-size red lid

JK171 Mini-size blue lid

JK172 Mini-size green lid

JK173 Mini-size gold lid

JK174 Mini-size silver lid



MD149 Mini-size silver lid with handles

Mini-Size Container Baskets

Mini Perforated Basket with Lid



JF159R

Mini-Size Tray for micro instruments fits into an Aesculap Mini container JK188 or JN188.

Length x Width x Height (in.)
Outside dimensions

JF159R 10½ x 5 x 2

Mini Perforated Basket with Lid, Silicone Cushioning Rack and Divider



JG380R

Mini-Size Tray contains three instrument racks fits into JN188 or JK188

Length x Width x Height (in.)
Outside dimensions

JG380R 10¾ x 3¾ x 1⅝

*All dimensions shown here are outside dimensions with the lid.

SterilContainer™

Wide-Body System

Wide-Body Containers – Length 23^{5⁄8} x Width 15^{3⁄4}

**Perforated Bottom
with Retention Plates**



JN817

4^{3⁄4}-inch height

Solid Bottoms



JK817



JN821

7-inch height



JK821

Wide-Body Container Lid



JE601

Wide-Body Basket



JF833R

Wide-Body Baskets

Length x Width x Height (in.)
Outside dimensions

JF833R 21 x 13^{1⁄4} x 3

*All dimensions shown here are outside dimensions with the lid.

JS Series SterilContainer™ S2

Container System Validated for Steam and Low Temperature

The Aesculap JS Series SterilContainer™ S2 Container System Product family is available for use with Low Temperature (STERRAD®, STERIS® V-PRO® and STERIZONE®), EtO and PreVac Steam sterilization modalities. See Aesculap Instructions for Use (IFU) www.aesculapusa.com.

Full-Size JS Series SterilContainer™ S2 Container System Bottoms – Length 23¼ x Width 11¼



Three-Quarter Size JS Series SterilContainer™ S2 Container System Bottoms – Length 18½ x Width 11¼




Half-Size JS Series SterilContainer™ S2 Container System Bottoms – Length 11¾ x Width 11¼



JS Series SterilContainer™ S2


Aluminum Lids

Full-Size, Three-Quarter and Half-Size JS Series SterilContainer™ S2 Container System Aluminum Lids



JS489

Full-Size Lid



JS789

Three-Quarter Size Lid



JS389

Half-Size Lid

SterilContainer™ S Mini-Size Container and Lid – Length 12¼ x Width 5½

Small heat and moisture-sensitive instruments can be secured in the JM188/JM174 Mini-Size container for low temperature sterilization and storage.



JM188

2⅝-inch height



JM174

Handle Option available as Ordering Code MD152

JM174 Lid WITHOUT handles (Shown)

MD152 Mini JM Container lid WITH handles (Picture Not Available)

Mini-Size Container Baskets

Mini Perforated Basket with Lid




JF159R

Mini-Size Tray for micro instruments fits into an Aesculap Mini container JK188 or JN188.

Length x Width x Height (in.)
Outside dimensions 10½ x 5 x 2

Perforated Silicone Holding Mats




JF944 – 9¾" x 4" Mini-size silicone mat

JF945 – 9¾" x 9¼" Half-size silicone mat

JF947 – 15¼" x 9¼" Three-Quarter size silicone mat

JF949 – 20¼" x 9¼" Full-Size silicone mat




SterilContainer™

System Processing Supplies

Tamper Evident Locks


Tamper Evident Locks Seal containers for safety and security. Indicator dots change color when exposed to sterilizing conditions. Quantities of 1000 each.

	Item No.	Color	Description	Qty.	PreVac Steam	PreVac IUSS	Gravity	EtO	Low Temp
	US900	Blue	No Indicator	1,000/pk	X	X	X	X	X
	US905	Green	Change ¹ Yellow to Orange	1,000/pk				X	
	US906	Orange	Change ¹ Blue to Brown	1,000/pk	X	X	X		
	US910	Pink	Change ¹ Magenta to Blue	1,000/pk					X ²
	US399	Yellow	Change ¹ Blue to Brown	1,000/pk		X			

1. After sterilization, the external indicator should change from the original indicator color to the recommended post sterilization indicator color. The recommended post sterilization indicator color may not be evenly shaded but should be shades of the proper processed color to indicate exposure to sterilant.

2. Locks must be stored in a controlled room temperature, away from alkaline chemicals, acids and sources of light.

Indicator Cards

	Item No.	Description	Qty.	PreVac Steam	PreVac IUSS	Gravity	EtO	Low Temp
	MD345	Communication Card No Indicator	250/pk	X	X	X	X	
	MD346	Change ^{1,2} Brown in Steam Change ^{1,2} Orange in EtO	250/pk	X	X	X	X	
	MD399	Change ¹ Brown in Steam	250/pk		X			
	MD876	Change ¹ Brown in Steam for Quarter and Mini Size	250/pk	X		X	X	
	MD334	Change ¹ Blue	250/pk					X
	MD335	Change ¹ Blue for JM188 with JM174	250/pk					X

1. After sterilization, the external indicator should change from the original indicator color to the recommended post sterilization indicator color. The recommended post sterilization indicator color may not be evenly shaded but should be shades of the proper processed color to indicate exposure to sterilant.

2. Filter contains a dual indicator dot, which changes to brown in steam, and to orange in EtO.

SterilContainer™

System Processing Supplies

Filters – Wide, Extra-Long, Full-Size, Three-Quarter and Half-Size Containers

	Item No.	Description	Qty.	PreVac Steam	PreVac IUSS	Gravity	EtO	Low Temp
	US751	Paper Filter with Indicator ¹	1,000/pk	X ¹	X ¹	X ¹	X ¹	
	US994	Paper Filter without Indicator	1,000/pk	X	X	X	X	
	MD344	Polypropylene Filter without Indicator	1,000/pk	X	X		X	X
	JK090	Metal Retention Plate Lid PTFE Reusable Filter	10/pk	X	X			
	JP050	PrimeLine™ PrimeLine Pro PTFE Reusable Filter	1/pk	X	X			

1. Filter contains a dual indicator dot, which changes from blue to brown in steam, and to orange in EtO.

Communication Cards, Biohazard Cards and Indicator Strips



US963 – Communicator Card
For use with US974 Communicator Tag.
Box of 500



US754 – Indicator Strip: (Class 1)
Single-use strip for use in steam and EtO sterilization. Box of 250



MD347 – Process Indicator Card
For use in steam or EtO sterilization. Box of 250 designed for Aesculap Generation 3 Containers, circa 1990's.



US991 – Biohazard Card
For use with soiled instruments.
Box of 250

*EDS Aicon version not approved for BioBin.


SterilContainer™

System Processing Supplies

Quarter-Size, Mini and Extra-Long Mini Size Containers

	Item No.	Description	Qty.	PreVac Steam	PreVac IUSS	Gravity	EtO	Low Temp
	US999	Paper Filter w/o Indicator	1,000/pk	X	X	X	X	
	MD355	Polypropylene Filter w/o Indicator	1,000/pk	X	X		X	X
	JK091	Metal Retention Plate Lid PTFE Reusable Filter	2/pk	X	X			

Filters – Generation 2 and Prior SterilContainer™ Systems



US992 — Single use
1/4-size
6 1/2" x 9" filter.
Box of 1000



JK092 — Reusable
6 1/2" x 9" filter.
Pack of 2



US756 — Single use
square 9" x 9" filter
with indicator dot.
Box of 1000



US998 — Single use
square 9" x 9" filter.
Box of 1000




JK089 — Reusable
9" x 9" filter.
Pack of 10

SterilContainer™

System Identification Labels and Tags

Identification Labels for Labeling Contents and Delivery Address



JG785U	Red
JG786U	Blue
JG787U	Green
JG788U	Gold
JG790U	Black
JG792U	Brown
JG793U	Purple
JG794U	Orange



ID Tag and Snap-on ID Tags US901, US917 and US945 are specific for the vapor hydrogen peroxide sterilization cycles.

Snap-on Identification Tags for Labeling Contents

Small Tags

2 3/8" x 1/2" (60 x 15 mm)

US935	Red
US936	Blue
US937	Green
US938	Gold
US939	Silver
US940	Black
US941	Purple
US942	Orange
US943	Bronze



Large Tags

2 3/8" x 1 1/8" (60 x 30 mm)

US965	Red
US966	Blue
US967	Green
US968	Gold
US969	Silver
US970	Black
US971	Purple
US972	Orange
US973	Bronze



Identification Labels for labeling contents for Vapor Hydrogen Peroxide Sterilization Systems.

US901	Magenta	US901
US917	Magenta	US917
US945	Magenta	US945



US901

ID Tag, 2 3/8" x 1/2"

US917

Small Snap-on ID Tag, 2 3/8" x 1/2"

US945

Large Snap-on ID Tag, 2 3/8" x 1 1/8"

To order this tag without clips add the suffix "200" behind the item #'s (example US935200 is the red small tag only)

Not recommended for use in vapor hydrogen peroxide sterilization due to color fading.

ID Tag and Snap-on ID Tags US901, US917 and US945 are specific for the vapor hydrogen peroxide sterilization cycles.

Instrument Organization System

Organization and protection of surgical instruments

Aesculap's Instrument Organization System is an advanced system of organizing, storing and protecting your instrument set during transportation and sterilization.

- Instruments organized with the Instrument Organization System are designed to guard against damage from instrument-on-instrument impact during transportation
- Instrument protection can help reduce the cost of repairs and set downtime
- Instrument protection can help increase instrument life and reduces replacement cost

In the O.R., a set organized with Aesculap's Instrument Organization System is designed to help cases flow more efficiently by allowing staff to:

- Quickly identify and select instruments
- Work directly from the organized set
- Replace used instruments into the slots making final instrument count easier


During sterilization, Aesculap's Instrument Organization System is designed for:

- Better sterilant flow between instruments
- Consistent sterilization




Bone tray accessory

Please contact your Aesculap Sales Representative or Customer Resource at (800) 282-9000 for a free consultation in organizing any of your surgical assets in the Aesculap SterilContainer™ System



Tympanoplasty set








Micro instrument tray

All set examples for illustrative purposes only.

Instrument Organization System


Silicone Instrument Racks (Mounting Type A)*

 JG304 – 4¾" x 1½"	 JG305 – 6¼" x 1½"	 JG306 – 9½" x 1½"	 JG399 – 6" x ½"
 JG308 – 4¾" x 1½"	 JG309 – 6¼" x 1½"	 JG310 – 9½" x 1½"	For basket cover  JG311 – 9½" x 1½"
 JG312 – 4¾" x 1½"	 JG313 – 6¼" x 1½"	 JG314 – 9½" x 1½"	
 JG316 – 4¾" x 1½"	 JG317 – 6¼" x 1½"	 JG318 – 9½" x 1½"	
 JG321 – 4⅞" x ¾"	 JG322 – 4⅞" x ¾"	 JG323 – 4⅞" x ¾"	 <i>Support for angled instruments</i> JG324 – 4⅞" x 1¾"
 JG361 – 6⅝" x 2"	PL962R  Trocar Storage Set Universal silicone storage racks for trocars		

*For all "Type A" mounting items, see Hardware and Tool Kit options on page 43.
JG303R option requires 7 mm Wrench/Socket and T15 sized screwdriver.
Hardware for mounting JG303R screws available as GB459R. See page 43.

Please contact your Aesculap Sales Representative or Customer Resource at (800) 282-9000 for a free consultation in organizing any of your surgical assets in the Aesculap SterilContainer™ System


Separating Racks (Mounting Type A)*




Small Silicone Holders (Mounting Type A)*



Metal Clamp with Silicone Inlay (Mounting Type A)*




Fixation Pins for Ring Handle Instruments



*For all "Type A" mounting items, see Hardware and Tool Kit options on page 43.

Instrument Organization System

Micro Instrument Organization System Instrument Racks




Large Silicone Holders (Mounting Type A)*



Mounting Instructions

Mounting A for instrument racks to be fixed with screws



MATERIALS REQUIRED:

- GB459R Toolkit as required
- Instrument holder such as JG341
- AESCULAP baskets (e.g. AESCULAP Aicon® JJ series or JF series)

Place the instrument holder in the required position in the AESCULAP basket, making sure that the holder's fixation lug fit onto the basket's lattice-work correctly. Push the screw through from below and fix it with the nut from above.

Hold the screw with your finger, put on the socket wrench with the other hand and tighten the nut.

*For all "Type A" mounting items, see Hardware and Tool Kit options on page 43.

Please contact your Aesculap Sales Representative or Customer Resource at (800) 282-9000 for a free consultation in organizing any of your surgical assets in the Aesculap SterilContainer™ System

Optics Holders (Mounting Type A)*

Eyepiece holder



JG367

Sleeve holder



JG368

Support for angled optics



JG369



Fixation Elements

Fixation pin for Mounting Type A
Can be applied using a flat head screwdriver



JG300 10 pins/pack

Socket Wrench



JJ589R

10 nuts and 10 bolts/pack



JJ588R

10 nuts and 10 bolts/pack



JG303R


Note: Bolt dimension is 7 mm and
Screw is T15

Mounting Hardware Tool Kit
for JG303R




GB459R

Tray Insert and Silicone Inserts (Mounting Type C)*



JF210R 17 1/4" x 10"

To be used with tray
inserts JF110F &
JF210F



JG370 9 1/2" x 1 1/8"

To be used with tray inserts
JF110F &
JF210F



JG371 9 1/2" x 1 1/8"

To be used with tray inserts
JF110F &
JF210F



JG372 9 1/2" x 1 1/8"

JF110R 10" x 7 3/4" - Designed for Half-Size JF Baskets

JF250R 12" x 10" - Designed for Three-Quarter Size JF Baskets


JF210R 17 1/4" x 10" - Designed for Full-Size JF Baskets

*For all "Type A" mounting items, see Hardware and Tool Kit options on page 43.

Instrument Organization System

Mounting Type A

For instrument racks with fixation pin:




Mounting Type A designed for easy top down assembly of the majority of SterilContainer™ Storage System silicone dividers

- ① Insert fixation pin JG300 in the rack. Install rack with fixation pin in the tray.
- ② Push the fixation pin down, then turn it by 45° using JG302 – Allen screwdriver.

Mounting Type B

For silicone holders with fixation clamps:




For use with JG330, JG331 and JG332 only

- ① Press fixation clamp JG301R into tray from below.
- ② Slide in silicone holder sideways.

Mounting Type C

For silicone racks in tray inserts:



For use when inserting JG310 into JF110R and JF210R

- ① Pull the silicone rack through the longitudinal slots in the tray.


**Mounting Type C can only be used with JF110R and JF210R baskets*

Please contact your Aesculap Sales Representative or Customer Resource at (800) 282-9000 for a free consultation in organizing any of your surgical assets in the Aesculap SterilContainer™ System

Instrument Holders



Item No.	Description
MB216	Holding Clamp I clamping capacity: 7-12 mm dia.
MB217	Holding Clamp II clamping capacity: 12-16 mm dia.
MB218	Holding Clamp III clamping capacity: 16-26 mm dia.
MB219	Holding Clamp IV clamping capacity: 16-36 mm dia.
MB220	Holding Clamp V clamping capacity: 36-45 mm dia.



Not for use in low temperature sterilization.

Holding Pins



*Holding Clamps and Holding Pins are compatible with both wire mesh perforated baskets.

Item No.	Description	Dimensions	
JF346	Holding Pin for Ring-Handled Instruments	1" dia. x 1 1/8" height	25 mm dia. x 29 mm height
US928	Holding Pin	1 1/4" height	28 mm height
US929	Holding Pin	1 3/4" height	45 mm height
US931	Holding Pin	2 1/4" height	56 mm height
US930	Holding Pin	3 1/2" height	90 mm height

Instrument Organization System


Micro Trays with Lid



JF149R – Mini-Half Size Basket
 $5\frac{3}{8}" \times 4\frac{3}{4}" \times 1\frac{7}{8}"$
 135 x 120 x 47 mm
Note: 2 Mini-Half Size Baskets fit into
 1 Mini Size SterilContainer™ System




**JF150R – Micro Instrument Basket
 with lid**
 $9\frac{1}{4}" \times 5\frac{3}{8}" \times 1\frac{5}{8}"$
 235 x 137 x 42 mm



Mini tray can be mounted on a
 perforated basket with
 fixation pin type
 A (JG300)

JF148R $5\frac{3}{4}" \times 3\frac{1}{4}" \times 1"$



Mini tray can be mounted on a
 perforated basket with
 fixation pin type
 A (JG300)

JF146R $3\frac{3}{8}" \times 1\frac{3}{4}" \times 1"$



Mini tray can be mounted on a
 perforated basket with
 fixation pin type
 A (JG300)

JF147R $14" \times 1\frac{3}{4}" \times 1"$

Microinstrument Trays




Components	Container Size		
	Full Size	Three-Quarter Size	Half Size
Bottoms	US437 – 21.25 x 10.00 x 1.50	US436 – 16.00 x 10.00 x 1.50	US435 – 10.00 x 9.88 x 1.50
Lids	US434 – 21.50 x 10.06 x 0.75	US433 – 16.50 x 10.88 x 0.75	US432 – 10.13 x 9.75 x 0.75

*For all "Type A" mounting items, see Hardware and Tool Kit options on page 43.

Please contact your Aesculap Sales Representative or Customer Resource at (800) 282-9000 for a free consultation in organizing any of your surgical assets in the Aesculap SterilContainer™ System

IOS Baskets – Some Assembly Required



JF449R

JF449R
Da Vinci Single Port (SP) Custom
Camera Container Insert



MD425

MD425 – Dual Scope Platform
(fits into Full-Size 5½" Container)
for both da Vinci Si series (without
cable connectors) and S series scopes
(with cable connectors)



MB071

MB071
Three-quarter sized custom basket for
microvascular instruments



MB072

MB072
Extra-Long custom basket for
Extra-Long robot instruments



MB073

MB073
Full-size custom basket for robot
instruments



MB074

MB074
Full-size
custom basket
for robot accessories



MB075

MB075
Full-size custom basket for
cystoscopes type 1




MB077

MB077
Full-size custom basket for hand
surgery instruments




MB078

MB078
Three-quarter
sized custom
basket for eye
surgery instruments



MB079

MB079
Full-size custom basket for cardiac
instruments




MB080

MB080
Three-quarter sized custom basket for
curettes



ME742

ME742
CASPAR Cervical Distractor (CCD)
Basket



US537

US537
Bipolar Forceps Basket
(bottom with lid)



MF320

MF320
Microdissection Instrument Basket




ME764

ME764
CASPAR Cervical Retractor (CCR)
Basket


Instrument Organization System

Set Examples


Transsphenoidal Hypophysectomy Set





Tympanoplasty Set




Laparoscopic Vascular Urology Set




Transnasal Pituitary Set




Micro Instrument Skull Base Set




Micro Cardio Vascular Instrument Set





Minimally Invasive Heart Valve Instrument Set



Trocar Set





CABG Micro Instrument Set




Instrument Organization System

Set Examples

Example Instrument Organization




Yellow Mats JF944, JF945, JF947, JF949
(fit into basket)






Example Instrument Organization



JG318



JG304



JG304




Call your local Aesculap Sales Representative for more set examples and onsite consultation


General Accessories

Endoscopy Storage Accessories


Endoscopy Storage Accessories: Scope Trays and Baskets




JF431R – Single Scope Tray
17⁷/₈" x 3³/₈" x 1⁵/₈"
454 x 84 x 41 mm



JF432R – Double Scope Tray
17⁷/₈" x 4¹/₁₆" x 1⁵/₈"
454 x 104 x 41 mm




JF433R – Double Scope Tray with Silicone Supports
17⁷/₈" x 4¹/₁₆" x 1⁵/₈"
454 x 104 x 41 mm



JF434R – Single Scope Tray with Silicone Rack
17⁷/₈" x 4¹/₁₆" x 1⁵/₈"
454 x 104 x 41 mm



JF436R – Single Scope Tray
(fits inside 1/4 size JN094)
10³/₄" x 3³/₈" x 1⁵/₈"
273 x 84 x 41 mm



JF435R – Scope Basket with Lid
(fits into Extra-Long Mini Container)
20" x 5¹/₄" x 1³/₄"
508 x 133 x 44 mm

All dimensions are length x width x height

Endorack Storage System for MIS Instruments



PL960R

PL960R
Endorack Full-Size,
Perforated basket and
silicone mesh mat, fits into
8" height, Full-Size SterilContainer™




PL961R

PL961R
Endorack Extra-Long,
Perforated basket and silicone mesh
mat, for instruments with shaft-
length up to 42 cm, fits into JN445
and also holds da Vinci SP Robotic
Arms



PL963R

PL963R
Endorack Full-Size,
fits into 8" height, Full-Size
SterilContainer



PL964R


PL964R
Endorack Extra-Long, fits into 8"
height, Extra-Long SterilContainer



PL960206


PL960206
Silicone Component
for PL960R, PL961R, PL963R and PL964R

Endoscopy Rack



MD379

Holds 14 Instruments



ME150


**ME150 – Full Size Multi-Keyhole
Bracket**
9 1/4" x 3/8" x 1 1/2"
236 x 9 x 38 mm

	Full-Size	Extra-Long
Rack Only	MD377	MD375
Rack with Basket	MD378 MD377 Rack, PL950215 Mat, JF222R Basket	MD379 MD375 Rack, JF933 Mat, JF232R Basket

General Accessories


Instrument Storage Racks

Racks



MD365

MD365 – Arthroscopy Rack
Holds 20 arthroscopy punches
Fits into a full-size, 8" deep
SterilContainer™
Also will fit with Solid Bottom
6" JK442, with JF222R
basket to hold additional
instrumentation associated
with punches
19 3/8" x 9 3/8" x 4 1/2"
490 x 238 x 114 mm




MD571

MD571 – Curette Tray
(bottom and lid)
Holds 14 curettes up to 6 1/4"
long



MD573

MD573 – Curette Tray
(bottom only)
Holds 14 curettes up to 6 1/4"
long



MD575

MD575 – Curette Tray
(lid only)
Holds 7 curettes up to 9" long



MD927

**MD927 – Osteotome Rack
for Three-Quarter Size
Container**
Holds 12 straight and curved
osteotomes (example: MB640
and MB650)
11 1/2" x 9 1/2" x 2 3/4"
295 x 240 x 70 mm



MD937

**MD937 – Osteotome Rack
for Half-Size Container**
Holds 12 straight and curved
osteotomes (example: MB640
and MB650)
10" x 9 7/8" x 2 3/4"
254 x 250 x 70 mm



MD928


MD928 – Cobb Curette Rack
Holds 10 Cobb curettes
(example: FK581 – FK599R)
9" x 2 1/2" x 14"
230 x 65 x 355 mm




MD929

**MD929 – Bruns Curette
Rack**
Holds 10 Bruns curettes
(example: MD473 – MD490)
9" x 2 1/2" x 11 1/2"
230 x 65 x 295 mm

Storage Frame for Bone Punches




JF120R holds ten instruments.
9½" x 8⅞" x 6⅜"



Instrument Set Example


JF120R can be used in the JN744, 8" three-quarter size container. It is designed to be used in a full-size container for longer instruments.

Lifting Platforms




Lifting Platform

Lifting Platform Material Code	Container Height	Handle Height*	Fits Container Codes
US094	4"	2½"	JK/JN 440
US095	5"	3½"	JK/JN 441
US096	6"	4¼"	JK/JN 442
US097	8"	6¼"	JK/JN 444
US098	10"	8½"	JK/JN 446



MD301

Fits into full-size 6" container JN/JK442 with handle folded 90°. Fits into full-size 10" container JN/JK446 with handles extended. 6¾" tall with handles extended vertically



JF101R

For added support for sets in a Full-Size, solid or perforated bottom container.

JF102R - Support for Half-Size

JF104R - Support for Three-Quarter Size



All dimensions are length x width x height

General Accessories

Basket Dividers

Basket Dividers



Item No.	Description	Outside Dimensions (length x width x height)	
JF378R	Limiting Bar with Securing Screws	6 1/2" x 1 1/16" x 1/2"	164 x 27 x 12 mm
JF381R	Limiting Bar with Securing Screws	3" x 1 1/16" x 1/2"	76 x 27 x 12 mm
JF921R	Straight divider for baskets	4 3/8" x 5/8" x 1 1/16"	110 x 15 x 28 mm
JF922R	Straight divider for baskets	6 3/8" x 5/8" x 1 1/16"	174 x 15 x 28 mm
JF923R	Straight divider for baskets	9 5/8" x 5/8" x 1 1/16"	244 x 15 x 28 mm
JG378	Basket Lid Instrument Support Holder	4 9/20" x 7/10" x 1 11/50"	113 x 18 x 21 mm
JG379	Basket Lid Instrument Support Holder	4 9/20" x 7/10" x 1 7/20"	113 x 18 x 31 mm

General Accessories

Micro Instruments Holders and Racks

Micro Instrument Holders and Racks



JF149R – Mini-Half Size Basket
5³/₈" x 4³/₄" x 1⁷/₈"
135 x 120 x 47 mm
Note: 2 Mini-Half Size Baskets fit into
1 Mini Size SterilContainer™ System



**JF150R – Micro Instrument Basket
with lid**
9¹/₄" x 5³/₈" x 1⁵/₈"
235 x 137 x 42 mm

Repair Tag




US990 – Repair Tag
1000/pack

General Accessories

Mats, Lubricants and Cleaners

Silicone Mesh Mats




JF938 – 9½" x 9⅞" (For use with half-size containers)

JF939 – 18⅞" x 9⅞"

JF940 – 15⅞" x 9⅞" (For use with three-quarter size containers)

JF941 – 21⅞" x 9⅞" (For use with full-size containers)

Silicone Cushioning Pads



MD894 Mini-size

Silicone Cushioning Pad

4⅜" x 9⅞"

110 x 250 mm

MD895 Quarter-size

Silicone Cushioning Pad

6½" x 10¼"

165 x 260 mm

MD896 Half-size

Silicone Cushioning Pad

9½" x 9⅞"

240 x 250 mm

MD897 Three-quarter-size

Silicone Cushioning Pad

15¾" x 9⅞"


400 x 250 mm

MD898 Full-size

Silicone Cushioning Pad

21¼" x 9⅞"

540 x 250 mm




JF933 Extra Long-size Silicone Cushioning Pad

24" x 9" x 1⅞"

610 x 230 x 30 mm

- Designed for safe storage of delicate instruments.
- Pad elevators designed to minimize potential moisture buildup.
- Perforated to help with air and steam circulation.





JF932 Silicone Mat, Blue – Full Size

18½" x 9" x 1" in

470 x 230 x 30 mm

All dimensions are length x width x height

Silicone Holding Mats with Perforations



- JF944 - 9¾" x 4" Mini-size silicone mats
- JF945 - 9¾" x 9¼" Half-size silicone mats
- JF947 - 15¼" x 9¼" Three-Quarter size silicone mat
- JF949 - 20¼" x 9¼" Full-Size silicone mat

Lubricants and Cleaners



JG598

JG598
Aesculap Instrument Oil for box locks, ratchets, etc.
 Physiologically inert acc. to DAB 7.



JG600

JG600
Aesculap Oil Spray for maintenance of
 instruments prior to sterilization. Physiologically
 inert acc. to DAB 7.
 Minimum order size 6, must be ordered in
 multiples of 6



JG601

JG601
Eloxal Cleaner, 100 ml for cleaning
 container surfaces.

General Accessories


Communicator Tags, Stringers and ID Brackets

Communicator Tag



US974 Communication Card Holder including Clip
 (Shown with US963 Communicator Card not included)
 Clip only - MD353
 US963 Dimensions: 3 1/2" x 2"

Stringers



All dimensions are from outside edge to outside edge

Notes: Approximately (+/-) 5 Ring Handled Instruments per each (1)" of Stringer length.

Open Stringers

US897	12" x 2 1/2"
US898	10" x 2 1/2"
US899	8" x 2 1/2"
US811	6" x 2 1/2"

Closed Stringers Push to Open


US890	12" x 2 1/2"
US889	10" x 2 1/2"
US888	8" x 2 1/2"
US926	6" x 2 1/2"
US814	4" x 2 1/2"
US909	12" x 5"
US925	10" x 5"
US904	8" x 5"

Bar Code ID Brackets for Baskets

JG616U and JG617U are Bar Code Identification brackets that can be attached to AESCULAP Aicon®, JJ Series, JF Series or perforated baskets to aid in securing bar code labels to the basket. The stainless steel brackets are secured to the baskets by bending the flexible tabs into the perforations.




58 x 23 mm (2 1/4" x 7/8")



89 x 30 mm (3 1/2" x 1 3/16")



96 x 34 mm (3 3/4" x 1 1/4")



105 x 54 mm (4 1/8" x 2 1/8")



Storage and Transport

Display Boards and Instrument Holders

Instrument Display Board



US872
Instrument Display Board only
34" x 34", Aluminum
864 x 864 mm
Picture shown with US863 Instrument Pegs, US862 Instrument Sleeves,
US866 Curette Holder, US865 Slotted Holder and various sized trays.
Each row (total 4) can accomadate (16) US863 Pegs

US867
Dust Cover

Instrument Holders and Board Trays



Item No.	Description
US865	Slotted Holder, 5" (127 mm)
US866	Curette Holder, 5" (127 mm) holds five curettes
US869	Instrument Board Tray, 11" (279 mm)
US870	Instrument Board Tray, 7" (178 mm)

Instrument Peg



US862
Instrument Sleeves
24/Pack

US863
Instrument Peg
6" (152 mm) (Qty = 1)

Miscellaneous

Metal Basins and Bowls

Metal Basins and Bowls



US060 - Metal Sponge Bowl
1½-quart capacity, 7" wide x 3" high



US061 - Metal Sponge Bowl
1-quart capacity, 6" wide x 2½" high



US062 - Metal Sponge Bowl
1⅞-quart capacity, 7½" wide x 3" high



US063 - Metal Iodine Cup
6 oz., 3½" wide x 2" high



US064 - Metal Graduated Measure
32 oz., 4½" wide x 5" high



US065 - Metal Iodine Cup
14 oz., 4½" wide x 2½" high



US066 - Metal Medicine Cup
2 oz., 2" wide x 2" high



US067 - Solution Basin
12 oz., 4½" wide x 2" high



JG504R - Kidney Basin
6¾", 250 ml, 1" high x 3½" wide



JG506R - Kidney Basin
10", 500 ml, 1" high x 5½" wide




JG508R - Kidney Basin
10¾", 750 ml, 1½" high x 6" wide

Miscellaneous

SterilContainer™ Replacement Parts

SterilContainer™ Replacement Parts



JK100 — Filter retainer
For use with lids and perforated container bottoms.



JK098 — Mini filter retainer



JP050 — Filter for PrimeLine™ lids and PrimeLine Pro



JP055 — Mounting handle for filter insert, PrimeLine lids



JF112210 — Replacement feet for perforated basket



JP001204 — Filter cover for PrimeLine and PrimeLine Pro



JP001202 — Outside cover for PrimeLine

Miscellaneous

Bio-Bin Containers

Bio-Bin Containers



The stainless steel *BIO-BIN* container is for transportation of soiled instruments and supplies within the hospital, designed to help prevent cross contamination. **Meets OSHA standards.**



JK262R



JK268R

JK262R
BIO-BIN Full-Size, Stainless Steel Bottom
23" x 11" x 6"
586 x 280 x 150 mm

JK268R
BIO-BIN Full-Size, Stainless Steel Lid
23" x 11" x 1 1/2"
586 x 280 x 36 mm



JK060R



JK062R

JK060R
BIO-BIN Half-Size, Stainless Steel Bottom
11 5/8" x 11" x 4"
295 x 280 x 102 mm

JK062R
BIO-BIN Half-Size, Stainless Steel Bottom
11 5/8" x 11" x 6"
295 x 280 x 150 mm



JK068R

JK068R
BIO-BIN Half-Size, Stainless Steel Lid
11 5/8" x 11" x 1 1/2"
295 x 280 x 36 mm

All dimensions are length x width x height

[illegible]

Sterile Container System Dimensions

TABLE OF DIMENSIONS FOR STERILE CONTAINERS SIZE FULL-SIZE
WITH AND WITHOUT ENHANCED DRYING SYSTEM (EDS)¹


Item no.	Net Weight (lbs)	Lid - outer dimensions (LxWxH)	Sterile container bottom - outer dimensions (LxWxH)	Sterile container bottom - inner dimensions [LxWxH] (loading height)	Sterile container with lid - outer dimensions (LxWxH)
Lid					
JJ410	4.114	23 ⁵ / ₈ x 11 ³ / ₄ x 1 ¹ / ₈ in 598 x 297 x 28 mm			
Sterile container bottom					
JJ110	3.523		23 ⁵ / ₈ x 11 ³ / ₄ x 3 ¹ / ₂ in	21 ⁵ / ₈ x 10 ¹ / ₄ x 2 ⁵ / ₈ in	23 ⁵ / ₈ x 11 ³ / ₄ x 4 in
JJ111	3.536		598 x 297 x 86 mm	548 x 258 x 64 mm	598 x 297 x 100 mm
JJ120	4.495		23 ⁵ / ₈ x 11 ³ / ₄ x 5 ³ / ₈ in	21 ⁵ / ₈ x 10 ¹ / ₄ x 4 ¹ / ₂ in	23 ⁵ / ₈ x 11 ³ / ₄ x 6 in
JJ121	4.508		598 x 297 x 136 mm	548 x 258 x 114 mm	598 x 297 x 150 mm
JJ130	5.467		23 ⁵ / ₈ x 11 ³ / ₄ x 7 ³ / ₈ in	21 ⁵ / ₈ x 10 ¹ / ₄ x 6 ¹ / ₂ in	23 ⁵ / ₈ x 11 ³ / ₄ x 8 in
JJ131	5.481		598 x 297 x 186 mm	548 x 258 x 164 mm	598 x 297 x 200 mm
JJ140	6.442		23 ⁵ / ₈ x 11 ³ / ₄ x 9 ³ / ₈ in	21 ⁵ / ₈ x 10 ¹ / ₄ x 8 ¹ / ₂ in	23 ⁵ / ₈ x 11 ³ / ₄ x 10 in
JJ141	6.455		598 x 297 x 236 mm	548 x 258 x 214 mm	598 x 297 x 250 mm

TABLE OF DIMENSIONS FOR STERILE CONTAINERS SIZE THREE-QUARTER
WITH AND WITHOUT ENHANCED DRYING SYSTEM (EDS)¹

Item no.	Net Weight (lbs)	Lid - outer dimensions (LxWxH)	Sterile container bottom - outer dimensions (LxWxH)	Sterile container bottom - inner dimensions [LxWxH] (loading height)	Sterile container with lid - outer dimensions (LxWxH)
Lid					
JJ430	1.626	18 ³ / ₈ x 11 ³ / ₄ x 1 ¹ / ₈ in 464 x 297 x 28 mm			
Sterile container bottom					
JJ310	2.818		18 ³ / ₈ x 11 ³ / ₄ x 3 ¹ / ₂ in	16 ³ / ₈ x 10 ¹ / ₄ x 2 ⁵ / ₈ in	18 ³ / ₈ x 11 ³ / ₄ x 4 in
JJ311	2.818		464 x 297 x 86 mm	414 x 258 x 64 mm	464 x 297 x 100 mm
JJ320	3.631		18 ³ / ₈ x 11 ³ / ₄ x 5 ³ / ₈ in	16 ³ / ₈ x 10 ¹ / ₄ x 4 ¹ / ₂ in	18 ³ / ₈ x 11 ³ / ₄ x 6 in
JJ321	3.631		464 x 297 x 136 mm	414 x 258 x 114 mm	464 x 297 x 150 mm
JJ330	4.445		18 ³ / ₈ x 11 ³ / ₄ x 7 ³ / ₈ in	16 ³ / ₈ x 10 ¹ / ₄ x 6 ¹ / ₂ in	18 ³ / ₈ x 11 ³ / ₄ x 8 in
JJ331	4.458		464 x 297 x 186 mm	414 x 258 x 164 mm	464 x 297 x 200 mm
JJ340	5.258		18 ³ / ₈ x 11 ³ / ₄ x 9 ³ / ₈ in	16 ³ / ₈ x 10 ¹ / ₄ x 8 ¹ / ₂ in	18 ³ / ₈ x 11 ³ / ₄ x 10 in
JJ341	5.271		464 x 297 x 236 mm	414 x 258 x 214 mm	464 x 297 x 250 mm

¹ Inches rounded up to closest 1/8 inch.

TABLE OF DIMENSIONS FOR STERILE CONTAINERS SIZE HALF-SIZE
WITH AND WITHOUT ENHANCED DRYING SYSTEM (EDS)¹

Item no.	Net Weight (lbs)	Lid - outer dimensions (LxWxH) 	Sterile container bottom - outer dimensions (LxWxH) 	Sterile container bottom - inner dimensions [LxWxH] (loading height) 	Sterile container with lid - outer dimensions (LxWxH) 
Lid					
JJ420	2.937	11 ⁷ / ₈ x 11 ³ / ₄ x 1 ¹ / ₈ in 300 x 297 x 28 mm			
Sterile container bottom					
JJ210	1.940		11 ⁷ / ₈ x 11 ³ / ₄ x 3 ¹ / ₂ in 300 x 297 x 86 mm	10 x 10 ¹ / ₈ x 2 ⁵ / ₈ in 250 x 258 x 64 mm	11 ⁷ / ₈ x 11 ³ / ₄ x 4 in 300 x 297 x 100 mm
JJ211	1.942				
JJ220	2.560		11 ⁷ / ₈ x 11 ³ / ₄ x 5 ³ / ₈ in 300 x 297 x 136 mm	10 x 10 ¹ / ₈ x 4 ¹ / ₂ in 250 x 258 x 114 mm	11 ⁷ / ₈ x 11 ³ / ₄ x 6 in 300 x 297 x 150 mm
JJ221	2.562				
JJ230	3.181		11 ⁷ / ₈ x 11 ³ / ₄ x 7 ³ / ₈ in 300 x 297 x 186 mm	10 x 10 ¹ / ₈ x 6 ¹ / ₂ in 250 x 258 x 164 mm	11 ⁷ / ₈ x 11 ³ / ₄ x 8 in 300 x 297 x 200 mm
JJ231	3.181				
JJ240	3.801		11 ⁷ / ₈ x 11 ³ / ₄ x 9 ³ / ₈ in 300 x 297 x 236 mm	10 x 10 ¹ / ₈ x 8 ¹ / ₂ in 250 x 258 x 214 mm	11 ⁷ / ₈ x 11 ³ / ₄ x 10 in 300 x 297 x 250 mm
JJ241	3.801				

The AESCULAP Aicon® Sterile Container System Instructions for Use are published on <http://www.aesculapusaifus.com/>

NOTES:

¹ Inches rounded up to closest 1/8 inch.

AESCULAP Aicon®

Sterile Container System Dimensions

TABLE OF DIMENSIONS FOR
FULL-SIZE BASKETS¹


Item no.	Basket lid outer dimensions (LxWxH)	Basket outer dimensions (LxWxH)	Basket inner dimensions (LxWxH)
			
Basket lid			
JJ167R	20 ⁷ / ₈ x 10 x ³ / ₄ in 528 x 253 x 18 mm		
Basket			
JJ161R		21 ¹ / ₄ x 10 x 1 ¹ / ₂ in 537 x 250 x 35 mm	21 ¹ / ₈ x 9 ⁷ / ₈ x 1 ³ / ₈ in 535 x 248 x 32 mm
JJ162R		21 ¹ / ₄ x 10 x 2 ¹ / ₄ in 537 x 250 x 55 mm	21 ¹ / ₈ x 9 ⁷ / ₈ x 2 ¹ / ₈ in 535 x 248 x 52 mm
JJ163R		21 ¹ / ₄ x 10 x 3 ¹ / ₄ in 537 x 250 x 80 mm	21 ¹ / ₈ x 9 ⁷ / ₈ x 3 ¹ / ₈ in 535 x 248 x 77 mm
JJ164R		21 ¹ / ₄ x 10 x 4 ¹ / ₄ in 537 x 250 x 105 mm	21 ¹ / ₈ x 9 ⁷ / ₈ x 4 ¹ / ₈ in 535 x 248 x 102 mm

TABLE OF DIMENSIONS FOR
THREE-QUARTER SIZE BASKETS¹

Item no.	Basket lid outer dimensions (LxWxH)	Basket outer dimensions (LxWxH)	Basket inner dimensions (LxWxH)
			
Basket lid			
JJ367R	15 ⁵ / ₈ x 10 x ³ / ₄ in 394 x 253 x 18 mm		
Basket			
JJ361R		15 ⁵ / ₈ x 10 x 1 ¹ / ₂ in 403 x 250 x 35 mm	15 ⁷ / ₈ x 10 x 1 ³ / ₈ in 401 x 248 x 32 mm
JJ362R		15 ⁵ / ₈ x 10 x 2 ¹ / ₄ in 403 x 250 x 55 mm	15 ⁷ / ₈ x 10 x 2 ¹ / ₈ in 401 x 248 x 52 mm
JJ363R		15 ⁵ / ₈ x 10 x 3 ¹ / ₄ in 403 x 250 x 80 mm	15 ⁷ / ₈ x 10 x 3 ¹ / ₈ in 401 x 248 x 77 mm
JJ364R		15 ⁵ / ₈ x 10 x 4 ¹ / ₄ in 403 x 250 x 105 mm	15 ⁷ / ₈ x 10 x 4 ¹ / ₈ in 401 x 248 x 102 mm

¹ Inches rounded up to closest 1/8 inch.

TABLE OF DIMENSIONS FOR
HALF-SIZE BASKETS¹

Item no.	Basket lid outer dimensions (L x W x H) 	Basket outer dimensions (L x W x H) 	Basket inner dimensions (L x W x H) 
Basket lid			
JJ267R	9 1/8 x 10 x 3/4 in 231 x 253 x 18 mm		
Basket			
JJ261R		9 1/2 x 10 x 1 1/2 in 240 x 250 x 35 mm	9 3/8 x 9 7/8 x 1 3/8 in 238 x 248 x 32 mm
JJ262R		9 1/2 x 10 x 2 1/4 in 240 x 250 x 55 mm	9 3/8 x 9 7/8 x 2 1/8 in 238 x 248 x 52 mm
JJ263R		9 1/2 x 10 x 3 1/4 in 240 x 250 x 80 mm	9 3/8 x 9 7/8 x 3 1/8 in 238 x 248 x 77 mm
JJ264R		9 1/2 x 10 x 4 1/4 in 240 x 250 x 105 mm	9 3/8 x 9 7/8 x 4 1/8 in 238 x 248 x 102 mm

NOTES:

¹ Inches rounded up to closest 1/8 inch.

SterilContainer™

System Product Dimensions

Full-Size Containers										
Length x Width x Height (Inches and mm)			Length x Width x Height (Inches and mm)							
JN440	Bottom Only		4.7 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		23¼ x 10¾ x 3½		592 x 274 x 90		23¼ x 11¼ x 4¼		592 x 285 x 108	
	Inside Dimension						21½ x 10⅞ x 2½		544 x 258 x 63	
JN441	Bottom Only		5.1 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		23¼ x 10¾ x 4¾		592 x 274 x 120		23¼ x 11¼ x 5½		592 x 285 x 138	
	Inside Dimension						21½ x 10⅞ x 3⅝		544 x 258 x 93	
JN442	Bottom Only		5.5 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		23¼ x 10¾ x 5¼		592 x 274 x 135		23¼ x 11¼ x 6		592 x 285 x 153	
	Inside Dimension						21½ x 10⅞ x 4¼		544 x 258 x 108	
JN444	Bottom Only		6.4 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		23¼ x 10¾ x 7⅜		592 x 274 x 187		23¼ x 11¼ x 8		592 x 285 x 205	
	Inside Dimension						21½ x 10⅞ x 6¼		544 x 258 x 160	
JN446	Bottom Only		7.4 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		23¼ x 10¾ x 9¾		592 x 274 x 247		23¼ x 11¼ x 10½		592 x 285 x 265	
	Inside Dimension						21½ x 10⅞ x 8⅝		544 x 258 x 220	
JK440	Bottom Only		4.2 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		23¼ x 10¾ x 3½		592 x 274 x 90		23¼ x 11¼ x 4¼		592 x 285 x 108	
	Inside Dimension						21½ x 10⅞ x 3		544 x 258 x 75	
JK441	Bottom Only		4.6 lbs		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		23¼ x 10¾ x 4¾		592 x 274 x 120		23¼ x 11¼ x 5½		592 x 285 x 138	
	Inside Dimension						21½ x 10⅞ x 4⅞		544 x 258 x 105	
JK442	Bottom Only		5.0 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		23¼ x 10¾ x 5¼		592 x 274 x 135		23¼ x 11¼ x 6		592 x 285 x 153	
	Inside Dimension						21½ x 10⅞ x 4¾		544 x 258 x 120	
JK444	Bottom Only		6.0 lbs		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		23¼ x 10¾ x 7⅜		592 x 274 x 187		23¼ x 11¼ x 8		592 x 285 x 205	
	Inside Dimension						21½ x 10⅞ x 6¾		544 x 258 x 172	
JK446	Bottom Only		6.9 lbs		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		23¼ x 10¾ x 9¾		592 x 274 x 247		23¼ x 11¼ x 10½		592 x 285 x 265	
	Inside Dimension						21½ x 10⅞ x 9⅞		544 x 258 x 232	

*Basket height includes feet

Perforated Baskets				
<i>Length x Width x Height (Inches)</i>			<i>Length x Width x Height (mm)</i>	
JF221R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
1.7 lbs.	21 ¹ / ₄ x 10 x 1 ¹ / ₂ *	21 ¹ / ₈ x 9 ⁷ / ₈ x ⁷ / ₈	540 x 253 x 36*	538 x 251 x 22
JF222R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
2.0 lbs.	21 ¹ / ₄ x 10 x 2 ¹ / ₄ *	21 ¹ / ₈ x 9 ⁷ / ₈ x 1 ⁵ / ₈	540 x 253 x 56*	538 x 251 x 42
JF223R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
2.4 lbs.	21 ¹ / ₄ x 10 x 3*	21 ¹ / ₈ x 9 ⁷ / ₈ x 2 ¹ / ₂	540 x 253 x 76*	538 x 251 x 62
JF224R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
2.8 lbs.	21 ¹ / ₄ x 10 x 4 ¹ / ₈ *	21 ¹ / ₈ x 9 ⁷ / ₈ x 3 ⁵ / ₈	540 x 253 x 106*	538 x 251 x 92

Basket Cover for Perforated Baskets		
<i>Length x Width x Height (Inches and mm)</i>		
JF227R	Outside (in.)	Outside (mm)
1.5 lbs.	21 ¹ / ₂ x 10 ¹ / ₁₆ x ³ / ₄	544 x 257 x 18

NOTES:

¹ Inches rounded up to closest 1/8 inch.

SterilContainer™

System Product Dimensions

Three-Quarter Size Containers				
<i>Length x Width x Height (Inches and mm)</i>			<i>Length x Width x Height (Inches and mm)</i>	
JN740	Bottom Only	3.8 lbs.	Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
<i>Outside Dimension</i>	18½ x 10¾ x 3½	470 x 274 x 90	18½ x 11¼ x 4¼	470 x 285 x 108
<i>Inside Dimension</i>			16¼ x 10⅛ x 2½	421 x 258 x 63
JN741	Bottom Only	4.3 lbs.	Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
<i>Outside Dimension</i>	18½ x 10¾ x 4¾	470 x 274 x 120	18½ x 11¼ x 5½	470 x 285 x 138
<i>Inside Dimension</i>			16¼ x 10⅛ x 3⅝	421 x 258 x 93
JN742	Bottom Only	4.5 lbs.	Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
<i>Outside Dimension</i>	18½ x 10¾ x 5¼	470 x 274 x 135	18½ x 11¼ x 6	470 x 285 x 153
<i>Inside Dimension</i>			16¼ x 10⅛ x 4¼	421 x 258 x 108
JN744	Bottom Only		Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
<i>Outside Dimension</i>	18½ x 10¾ x 7⅝	470 x 274 x 187	18½ x 11¼ x 8	470 x 285 x 205
<i>Inside Dimension</i>			16¼ x 10⅛ x 6¼	421 x 258 x 172
JK740	Bottom Only	3.6 lbs.	Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
<i>Outside Dimension</i>	18½ x 10¾ x 3½	470 x 274 x 90	18½ x 11¼ x 4¼	470 x 285 x 108
<i>Inside Dimension</i>			16¼ x 10⅛ x 3	421 x 258 x 75
JK741	Bottom Only	4.0 lbs.	Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
<i>Outside Dimension</i>	18½ x 10¾ x 4¾	470 x 274 x 120	18½ x 11¼ x 5½	470 x 285 x 138
<i>Inside Dimension</i>			16¼ x 10⅛ x 4⅛	421 x 258 x 105
JK742	Bottom Only	4.2 lbs.	Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
<i>Outside Dimension</i>	18½ x 10¾ x 5¼	470 x 274 x 135	18½ x 11¼ x 6	470 x 285 x 153
<i>Inside Dimension</i>			16¼ x 10⅛ x 4¾	421 x 258 x 120
JK744	Bottom Only		Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
<i>Outside Dimension</i>	18½ x 10¾ x 7⅝	470 x 274 x 187	18½ x 11¼ x 8	470 x 285 x 205
<i>Inside Dimension</i>			16¼ x 10⅛ x 6¾	421 x 258 x 172

*Basket height includes feet

Perforated Baskets				
<i>Length x Width x Height (Inches)</i>			<i>Length x Width x Height (mm)</i>	
JF251R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
1.4 lbs.	16 x 10 x 1½*	15 ⁷ / ₈ x 9 ⁷ / ₈ x ⁷ / ₈	406 x 253 x 36*	404 x 251 x 22
JF252R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
1.6 lbs.	16 x 10 x 2¼*	15 ⁷ / ₈ x 9 ⁷ / ₈ x 1 ⁵ / ₈	406 x 253 x 56*	404 x 251 x 42
JF253R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
2.0 lbs.	16 x 10 x 3*	15 ⁷ / ₈ x 9 ⁷ / ₈ x 2½	406 x 253 x 76*	404 x 251 x 62
JF254R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
2.3 lbs.	16 x 10 x 4½*	15 ⁷ / ₈ x 9 ⁷ / ₈ x 3 ⁵ / ₈	406 x 253 x 106*	404 x 251 x 92

Basket Cover for Perforated Baskets		
<i>Length x Width x Height (Inches and mm)</i>		
JF257R	Outside (in.)	Outside (mm)
1.2 lbs.	16½ x 10 ⁷ / ₈ x ¾	419 x 276 x 19

NOTES:

* Inches rounded up to closest 1/8 inch.

SterilContainer™

System Product Dimensions

Half-Size Containers										
Length x Width x Height (Inches and mm)			Length x Width x Height (Inches and mm)							
JN340	Bottom Only		3.0 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		11 ³ / ₄ x 10 ³ / ₄ x 3 ¹ / ₂		300 x 274 x 90		11 ³ / ₄ x 11 ¹ / ₄ x 4 ¹ / ₄		300 x 285 x 108	
	Inside Dimension						10 x 10 ¹ / ₈ x 2 ¹ / ₂		254 x 258 x 63	
JN341	Bottom Only		3.3 lbs		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		11 ³ / ₄ x 10 ³ / ₄ x 4 ³ / ₄		300 x 274 x 120		11 ³ / ₄ x 11 ¹ / ₄ x 5 ¹ / ₂		300 x 285 x 138	
	Inside Dimension						10 x 10 ¹ / ₈ x 3 ⁵ / ₈		254 x 258 x 93	
JN342	Bottom Only		3.5 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		11 ³ / ₄ x 10 ³ / ₄ x 5 ¹ / ₄		300 x 274 x 135		11 ³ / ₄ x 11 ¹ / ₄ x 6		300 x 285 x 153	
	Inside Dimension						10 x 10 ¹ / ₈ x 4 ¹ / ₄		254 x 258 x 108	
JN344	Bottom Only		4.0 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		11 ³ / ₄ x 10 ³ / ₄ x 7 ³ / ₈		300 x 274 x 187		11 ³ / ₄ x 11 ¹ / ₄ x 8		300 x 285 x 205	
	Inside Dimension						10 x 10 ¹ / ₈ x 6 ¹ / ₄		254 x 258 x 160	
JN346	Bottom Only		4.7 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		11 ³ / ₄ x 10 ³ / ₄ x 9 ³ / ₄		300 x 274 x 247		11 ³ / ₄ x 11 ¹ / ₄ x 10 ¹ / ₂		300 x 285 x 265	
	Inside Dimension						10 x 10 ¹ / ₈ x 8 ⁵ / ₈		254 x 258 x 220	
JK340	Bottom Only		2.6 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		11 ³ / ₄ x 10 ³ / ₄ x 3 ¹ / ₂		300 x 274 x 90		11 ³ / ₄ x 10 ³ / ₄ x 4 ¹ / ₄		300 x 285 x 108	
	Inside Dimension						10 x 10 ¹ / ₈ x 3		254 x 258 x 75	
JK341	Bottom Only		3.1 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		11 ³ / ₄ x 10 ³ / ₄ x 4 ³ / ₄		300 x 274 x 120		11 ³ / ₄ x 10 ³ / ₄ x 5 ¹ / ₂		300 x 285 x 138	
	Inside Dimension						10 x 10 ¹ / ₈ x 4 ¹ / ₄		254 x 258 x 108	
JK342	Bottom Only		3.2 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		11 ³ / ₄ x 10 ³ / ₄ x 5 ¹ / ₄		300 x 274 x 135		11 ³ / ₄ x 10 ³ / ₄ x 6		300 x 285 x 153	
	Inside Dimension						10 x 10 ¹ / ₈ x 4 ³ / ₄		254 x 258 x 120	
JK344	Bottom Only		3.9 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		11 ³ / ₄ x 10 ³ / ₄ x 7 ³ / ₈		300 x 274 x 187		11 ³ / ₄ x 10 ³ / ₄ x 8		300 x 285 x 205	
	Inside Dimension						10 x 10 ¹ / ₈ x 6 ³ / ₄		254 x 258 x 172	
JK346	Bottom Only		4.5 lbs.		Bottom and Aluminum Lid					
	Inches		mm		Inches		mm			
	Outside Dimension		11 ³ / ₄ x 10 ³ / ₄ x 9 ³ / ₄		300 x 274 x 247		11 ³ / ₄ x 10 ³ / ₄ x 10 ¹ / ₂		300 x 285 x 265	
	Inside Dimension						10 x 10 ¹ / ₈ x 9 ¹ / ₈		254 x 258 x 232	

* Basket height includes feet

Perforated Baskets				
<i>Length x Width x Height (Inches)</i>			<i>Length x Width x Height (mm)</i>	
JF111R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
1.1 lbs.	10 x 9 ⁷ / ₈ x 1 ¹ / ₂ *	9 ¹ / ₂ x 9 ⁵ / ₈ x ⁷ / ₈	253 x 243 x 36*	251 x 241 x 22
JF112R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
1.3 lbs.	10 x 9 ⁷ / ₈ x 2 ¹ / ₄ *	9 ¹ / ₂ x 9 ⁵ / ₈ x 1 ⁵ / ₈	253 x 243 x 56*	251 x 241 x 42
JF113R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
1.5 lbs.	10 x 9 ⁷ / ₈ x 3*	9 ¹ / ₂ x 9 ⁵ / ₈ x 2 ¹ / ₂	253 x 243 x 76*	251 x 241 x 62
JF114R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
1.8 lbs.	10 x 9 ⁷ / ₈ x 4 ¹ / ₈ *	9 ¹ / ₂ x 9 ⁵ / ₈ x 3 ⁵ / ₈	253 x 243 x 106*	251 x 241 x 92

Basket Cover for Perforated Baskets		
<i>Length x Width x Height (Inches and mm)</i>		
JF117R	Outside (in.)	Outside (mm)
.8 lbs.	10 ¹ / ₈ x 9 ³ / ₄ x ³ / ₄	257 x 247 x 18

NOTES:

¹ Inches rounded up to closest 1/8 inch.

SterilContainer™

System Product Dimensions

Quarter-Size Containers for Dental

Filter Requirements: Perforated Bottoms with lids require two US992 filters. (6.5" x 9")

<i>Length x Width x Height (Inches)</i>			<i>Length x Width x Height (mm)</i>	
JN092	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
1.8 lbs.	12 ³ / ₄ x 7 ¹ / ₂ x 1 ⁷ / ₈	11 x 7 x ³ / ₄	313 x 190 x 46	277 x 177 x 20
JN094	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
1.8 lbs.	12 ³ / ₄ x 7 ¹ / ₂ x 2 ¹ / ₂	11 x 7 x 1 ¹ / ₂	313 x 190 x 65	277 x 177 x 39
JN095	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
	12 ³ / ₄ x 7 ¹ / ₂ x 3 ⁵ / ₈	11 x 7 x 2 ¹ / ₂	313 x 190 x 92	277 x 177 x 66
JN096	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
1.9 lbs.	12 ³ / ₄ x 7 ¹ / ₂ x 5 ¹ / ₈	11 x 7 x 4	313 x 190 x 130	277 x 177 x 104

Quarter-Size Lid Only

<i>Length x Width x Height (Inches and mm)</i>		
JN091	Outside (in.)	Outside (mm)
	12 ¹ / ₂ x 7 ¹ / ₂ x 1	319 x 190 x 25

Quarter-Size Container Baskets

Length x Width x Height (Inches)			Length x Width x Height (mm)	
JF264R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
1.2 lbs.	10 ³ / ₄ x 7 x 1 ⁵ / ₈	10 ¹ / ₈ x 6 ⁵ / ₈ x 1 ³ / ₈	274 x 177 x 42	258 x 170 x 34
Quarter Size Perforated Basket with Lid and Silicone Cushioning Rack				
JF282R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
1.5 lbs.	10 ⁵ / ₈ x 6 ³ / ₄ x 3 ⁵ / ₈	10 ¹ / ₂ x 6 ⁵ / ₈ x 3 ¹ / ₂	270 x 171 x 93	267 x 169 x 91
JF289R	Outside (in.)	Outside (mm)		
.6 lbs.	10 ³ / ₄ x 7 x ⁵ / ₈	273 x 177 x 16		
Quarter Size Lid for tray JF282R				
JG380R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
	10 ³ / ₄ x 3 ⁵ / ₁₆ x 1 ³ / ₁₆	10 ¹ / ₈ x 3 x 1 ¹ / ₁₆	273 x 84 x 30	257 x 76 x 27
Mini-Size perforated basket with lid and three instrument racks, two side by side wil fit into JN094				
JG381R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
	10 ³ / ₄ x 7 x 1 ³ / ₁₆	10 ¹ / ₈ x 6 ⁵ / ₈ x 1 ¹ / ₁₆	273 x 177 x 30	257 x 168 x 27
Quarter-Size perforated basket with lid				
JG384R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
	10 ³ / ₄ x 7 x 1 ⁵ / ₈	10 ¹ / ₈ x 6 ⁵ / ₈ x 1 ¹ / ₂	273 x 177 x 41	257 x 168 x 38
Quarter-Size perforated basket with lid and three instrument racks, one JG381R and one JG384R can fit into JN095				
MD374	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
	10 ³ / ₄ x 7 x 1 ⁵ / ₈	10 ¹ / ₈ x 6 ⁵ / ₈ x 1 ³ / ₈	274 x 177 x 42	258 x 170 x 34

Quarter Size Perforated Basket with Lid and Silicone Cushioning Rack

*Basket height includes feet

Mini-Size Containers				
<i>Length x Width x Height (Inches and mm)</i>			<i>Length x Width x Height (Inches and mm)</i>	
JN187	Bottom Only	.5 lbs.	Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
<i>Outside Dimension</i>	12 ¹ / ₄ x 5 ¹ / ₈ x 1 ¹ / ₈	310 x 132 x 30	12 ¹ / ₄ x 5 ¹ / ₂ x 1 ¹ / ₂	310 x 140 x 40
<i>Inside Dimension</i>			10 ³ / ₄ x 5 x 1 ¹ / ₄	274 x 126 x 33
JN188	Bottom Only	.6 lbs.	Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
<i>Outside Dimension</i>	12 ¹ / ₄ x 5 ¹ / ₈ x 2 ¹ / ₄	310 x 132 x 57	12 ¹ / ₄ x 5 ¹ / ₂ x 2 ⁵ / ₈	310 x 140 x 67
<i>Inside Dimension</i>			10 ³ / ₄ x 5 x 2	274 x 126 x 50
JK187	Bottom Only	.4 lbs.	Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
<i>Outside Dimension</i>	12 ¹ / ₄ x 5 ¹ / ₈ x 1 ¹ / ₈	310 x 132 x 30	12 ¹ / ₄ x 5 ¹ / ₂ x 1 ¹ / ₂	310 x 140 x 40
<i>Inside Dimension</i>			10 ³ / ₄ x 5 x 1	274 x 126 x 25
JK188	Bottom Only	.6 lbs.	Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
<i>Outside Dimension</i>	12 ¹ / ₄ x 5 ¹ / ₈ x 2 ¹ / ₄	310 x 132 x 57	12 ¹ / ₄ x 5 ¹ / ₂ x 2 ⁵ / ₈	310 x 140 x 67
<i>Inside Dimension</i>			10 ³ / ₄ x 5 x 2	274 x 126 x 52

Mini-Size Container Baskets				
<i>Length x Width x Height (Inches)</i>			<i>Length x Width x Height (mm)</i>	
JF159R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
	10 ¹ / ₂ x 5 x 2	10 ¹ / ₄ x 4 ⁵ / ₈ x 1 ⁵ / ₈	267 x 125 x 49	261 x 118 x 41
JG380R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
	10 ³ / ₄ x 3 ³ / ₈ x 1 ⁵ / ₈	10 ¹ / ₂ x 3 ¹ / ₈ x 1 ³ / ₈	273 x 84 x 41	267 x 77 x 33

NOTES:

¹ Inches rounded up to closest 1/8 inch.

SterilContainer™

System Product Dimensions

Extra-Long Containers				
Length x Width x Height (Inches and mm)			Length x Width x Height (Inches and mm)	
JN445	Bottom Only	7.3 lbs.	Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
Outside Dimension	28 x 10 ³ / ₄ x 7 ³ / ₈	711 x 274 x 187	28 x 11 x 8	711 x 280 x 205
Inside Dimension			26 ¹ / ₄ x 10 ¹ / ₄ x 6 ¹ / ₄	667 x 261 x 160

Extra-Long Perforated Basket				
Length x Width x Height (Inches)			Length x Width x Height (mm)	
JF232R	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
	26 x 10 x 2 ¹ / ₄	25 ⁷ / ₈ x 9 ⁷ / ₈ x 1 ⁵ / ₈	662 x 254 x 56	658 x 251 x 42

Wide-Body containers				
Length x Width x Height (Inches)			Length x Width x Height (mm)	
JN817	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
	23 ⁵ / ₈ x 15 ³ / ₄ x 4 ³ / ₄	21 ¹ / ₂ x 13 ³ / ₄ x 3 ⁵ / ₈	600 x 400 x 120	544 x 350 x 91
JN821	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
	23 ⁵ / ₈ x 15 ³ / ₄ x 7	21 ¹ / ₂ x 13 ³ / ₄ x 6	600 x 400 x 180	544 x 350 x 151
JK817	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
	23 ⁵ / ₈ x 15 ³ / ₄ x 4 ³ / ₄	21 ¹ / ₂ x 13 ³ / ₄ x 4	600 x 400 x 120	544 x 350 x 103
JK821	Outside (in.)	Inside (in.)	Outside (mm)	Inside (mm)
	23 ⁵ / ₈ x 15 ³ / ₄ x 7	21 ¹ / ₂ x 13 ³ / ₄ x 6 ¹ / ₂	600 x 400 x 180	544 x 350 x 163

Wide-Body Baskets				
Length x Width x Height (Inches and mm)				
JE167R	Outside (in.)	Outside (mm)		
3.1 lbs.	21 x 13 ¹ / ₄ x 2	533 x 337 x 51		
JE168R	Outside (in.)	Outside (mm)		
3.7 lbs.	21 x 13 ¹ / ₄ x 3	533 x 337 x 76		
JE169R	Outside (in.)	Outside (mm)		
4.0 lbs.	21 x 13 ¹ / ₄ x 4	533 x 337 x 102		

*Basket height includes feet

SterilContainer™ S				
Length x Width x Height (Inches and mm)			Length x Width x Height (Inches and mm)	
JM440	Bottom Only		Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
Outside Dimension	23 ¹ / ₄ x 10 ³ / ₄ x 3 ¹ / ₂	592 x 274 x 90	23 ¹ / ₄ x 11 ¹ / ₄ x 4 ¹ / ₄	592 x 285 x 108
Inside Dimension			21 ¹ / ₂ x 10 ¹ / ₈ x 2 ¹ / ₂	544 x 258 x 63
JM441	Bottom Only		Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
Outside Dimension	23 ¹ / ₄ x 10 ³ / ₄ x 4 ³ / ₄	592 x 274 x 120	23 ¹ / ₄ x 11 ¹ / ₄ x 5 ¹ / ₂	592 x 285 x 138
Inside Dimension			21 ¹ / ₂ x 10 ¹ / ₈ x 3 ⁵ / ₈	544 x 258 x 93
JM442	Bottom Only		Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
Outside Dimension	23 ¹ / ₄ x 10 ³ / ₄ x 5 ¹ / ₄	592 x 274 x 135	23 ¹ / ₄ x 11 ¹ / ₄ x 6	592 x 285 x 153
Inside Dimension			21 ¹ / ₂ x 10 ¹ / ₈ x 4 ¹ / ₄	544 x 258 x 108
JM740	Bottom Only		Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
Outside Dimension	18 ¹ / ₂ x 10 ³ / ₄ x 3 ¹ / ₂	470 x 274 x 90	18 ¹ / ₂ x 11 ¹ / ₄ x 4 ¹ / ₄	470 x 285 x 108
Inside Dimension			16 ¹ / ₄ x 10 ¹ / ₈ x 2 ¹ / ₂	421 x 258 x 63
JM741	Bottom Only		Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
Outside Dimension	18 ¹ / ₂ x 10 ³ / ₄ x 4 ³ / ₄	470 x 274 x 120	18 ¹ / ₂ x 11 ¹ / ₄ x 5 ¹ / ₂	470 x 285 x 138
Inside Dimension			16 ¹ / ₄ x 10 ¹ / ₈ x 3 ⁵ / ₈	421 x 258 x 93
JM742	Bottom Only		Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
Outside Dimension	18 ¹ / ₂ x 10 ³ / ₄ x 5 ¹ / ₄	470 x 274 x 135	18 ¹ / ₂ x 11 ¹ / ₄ x 6	470 x 285 x 153
Inside Dimension			16 ¹ / ₄ x 10 ¹ / ₈ x 4 ¹ / ₄	421 x 258 x 108
JM340	Bottom Only		Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
Outside Dimension	11 ³ / ₄ x 10 ³ / ₄ x 3 ¹ / ₂	300 x 274 x 90	11 ³ / ₄ x 11 ¹ / ₄ x 4 ¹ / ₄	300 x 285 x 108
Inside Dimension			10 x 10 ¹ / ₈ x 2 ¹ / ₂	254 x 258 x 63
JM341	Bottom Only		Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
Outside Dimension	11 ³ / ₄ x 10 ³ / ₄ x 4 ³ / ₄	300 x 274 x 120	11 ³ / ₄ x 11 ¹ / ₄ x 5 ¹ / ₂	300 x 285 x 138
Inside Dimension			10 x 10 ¹ / ₈ x 3 ⁵ / ₈	254 x 258 x 93
JM342	Bottom Only		Bottom and Aluminum Lid	
	Inches	mm	Inches	mm
Outside Dimension	11 ³ / ₄ x 10 ³ / ₄ x 5 ¹ / ₄	300 x 274 x 135	11 ³ / ₄ x 11 ¹ / ₄ x 6	300 x 285 x 153
Inside Dimension			10 x 10 ¹ / ₈ x 4 ¹ / ₄	254 x 258 x 108

Through collaborative excellence we will improve the quality of a patient's life and meet the needs of the changing healthcare environment.



Instructions For Use are located at
www.aesculapusaifus.com



*Instructional videos on processing the
Aesculap SterilContainer™ System are located at*
www.youtube.com/user/aesculapusa

Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogs and price lists and to take legal actions.

©2024 AESCULAP, INC. ALL RIGHTS RESERVED. PRINTED IN THE USA.
Aesculap is an equal opportunity employer

Aesculap, Inc. | 3773 Corporate Parkway | Center Valley, PA | 18034
800-282-9000 | aesculapusa.com

Aesculap, Inc. – a member of the B. Braun Group of Companies

DOC132 Rev L_ 08/24