



Specialty Non-Stick Disposable Bipolar Forceps

Irrigating and Non-Irrigating

ENVISION THE EXTRAORDINARY™



AESCULAP®

Specialty Non-Stick Disposable Bipolar Forceps

Irrigating and Non-Irrigating

Features and benefits:

- Reliable single-use design helps ensure consistent performance and minimizes operating room disruptions
- Ergonomic handle designed for balance and precision
- Non-stick tips help reduce tissue adhesion
- Easy setup with integrated 12-foot cord (non-irrigating only)

**49%
HIGHER**
output than
competitive dual
irrigation¹



Irrigating

Item No.	Length	Tip Size
US981SU	7"	0.5 mm
US982SU	7"	1.0 mm
US983SU	7"	1.5 mm
US984SU	8"	0.5 mm
US985SU	8"	1.0 mm
US986SU	8"	1.5 mm
US987SU	9"	0.5 mm
US988SU	9"	1.0 mm
US989SU	9"	1.5 mm



Non-Irrigating

Item No.	Length	Tip Size
US801SU	7"	0.5 mm
US802SU	7"	1.0 mm
US803SU	7"	1.5 mm
US804SU	8"	0.5 mm
US805SU	8"	1.0 mm
US806SU	8"	1.5 mm
US807SU	9"	0.5 mm
US808SU	9"	1.0 mm
US809SU	9"	1.5 mm

Indications for Use:

Bipolar coagulation forceps are used for hemostatic coagulation in surgical procedures. See Instructions for Use for additional information, including warnings and precautions.

Aesculap, Inc. - a member of the B. Braun Group of Companies
Center Valley, PA | 800-282-9000 | aesculapusa.com