


DS Ligation Clip System

For Laparoscopic Appendectomy



Aesculap Surgical Technologies – Endoscopy

DS Ligation Clip System





Features & Benefits

Cost effective

- Up to 75% per procedure cost savings compared to endomechanical staplers

Diamond-shaped clip texture

- Maximizes contact surface on tissue to further enhance clip security, while remaining atraumatic

"Double Clip" design

- Acts like a fenestrated grasper for a firm hold on tissue and to stabilize the clip against axial (side-to-side) dislodgement


Tooth & notch

- Stabilize the clip against radial (pull off) dislodgement

Distal to proximal clip closure

- Prevents pushing out tissue during closure





PL809R
Endoscopic Clip Applier, X-large,
diam. 12 mm, length 310 mm




PL475SU
Titanium ligating clip, X-large
12 cartridges, 4 clips per cartridge

DS Ligation Clip System

Delivering cost savings without the compromise

The DS Ligation Clip System provides a cost-effective solution, without compromising on security, for laparoscopic appendectomies. This is made possible by an innovative extra large clip designed especially with closure of the appendix and mesoappendix in mind. The DS Ligation Clip System provides surgeons with **double** the clip for **savings, security** and **simplicity**.



DS Ligation Clips applied to appendix and mesoappendix

Indication for Use:

Intended for use in endoscopic surgery for ligating and marking vessels and tubular structures whenever ligating clips are used/indicated.

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogs and price lists and to take legal actions.

©2011 AESCULAP. ALL RIGHTS RESERVED. PRINTED IN THE USA.
Aesculap is an equal opportunity employer

Aesculap, Inc. | 3773 Corporate Parkway | Center Valley, PA | 18034
Phone 800-282-9000 | Fax 610-791-6886 | www.aesculapusa.com

Aesculap, Inc. - a B. Braun company

DOC1012 1M 2/11