

Laparoscopic Needleholders

Suturing Made Simple



Aesculap Surgical Technologies — Endoscopy

AESCULAP®


Suturing Made Simple



Balance





Control



Power

Aesculap's Axial Needle Holders – Suturing Made Simple

- **Balance** - revolutionary handle design provides superior balance
- **Control** - when closed, the handle is circular in the hand allowing precise movements and providing excellent ergonomics
- **Power** - force transmission to the jaws which grasps and holds needles securely





















Specialty "Super Righting"

The super righting needleholder saves time by automatically positioning the needle upright, while allowing surgeons steady control of the needle's entry.



Unparalleled Grasping

Even the finest needles and thinnest suture won't slip from these powerful, precision-engineered jaws.

Intracorporeal		31 cm	37 cm	42 cm	
	Durogrip TC Needle Holder, straight, 5 mm	PL407R	PL404R	PL416R	
	Durogrip TC Needle Holder, curved left, 5 mm	PL408R	PL405R	PL417R	
	Durogrip TC Needle Holder, curved right, 5 mm	PL413R	PL415R		
	"The Viper" Durogrip TC Suture Assistant /Needle Holder, 2 x 3 teeth, 5 mm	PL409R			
	Super Righting Needle Holder, 5 mm	PL414R			
	Durogrip TC Needle Holder, Fenestrated, 5 mm	PL410R			
	Durogrip TC Needle Holder, Crimped/Fenestrated, 5 mm	PL411R			
Extracorporeal					
	CROCE/OLMI Dissecting and Suture Scissors with Knot Pusher, 5 mm	PO022R			
	*Knot Pushing/Knot Tightening Forceps, 5 mm	PM654R PM973R PO958R	Jaw Insert Outer Tube Handle		
		PL651R			
		PL432R			
Reusable Scissors					
	PO885SU	Disposable Slightly Curved Metzenbaum Scissors Jaw Insert, 5 mm, 31 cm, 10 per box			
	PO886	Disposable Hook Scissors Jaw Insert, 5 mm, 31 cm, 10 per box			
	PO887	Disposable Mini-Metz Scissors, 5 mm, 31 cm, 10 per box			
	PO888	Disposable Metzenbaum Scissors Jaw Insert, 5 mm, 31 cm, 10 per box			
	PO889	Disposable Metzenbaum Scissors Jaw Insert, 5 mm, 42 cm, 10 per box			
	PO958R	Quick-Snap™ Handle		PM950R	Insulated Metal Handle

* Instruments sold by component

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogs and price lists and to take legal actions.

©2011 AESCULAP. ALL RIGHTS RESERVED. PRINTED IN THE USA.
Aesculap is an equal opportunity employer

Aesculap, Inc. | 3773 Corporate Parkway | Center Valley, PA | 18034
Phone 800-282-9000 | Fax 610-791-6886 | www.aesculapusa.com

Aesculap, Inc. – a B. Braun company

DOC460 Rev. C 5M 5/11