


ELAN[®] 4 Tool Catalog

Precise, Powerful, Intuitive.



Aesculap Power Systems

AESCULAP

ELAN[®] 4

Tool Catalog

Table of Contents

› Sagittal Saw Blades (Oscillating Saw Blades)	3-7
› Reciprocating Saw Blades	8-12
› Transversal Saw Blades	13-14
› ELAN 4 Burrs	15-23


All measurements and technical data listed in this brochure are approximate, as they are subject to normal manufacturing tolerances and technical modifications.

Indications for USE:

The ELAN 4 Tools are intended for high speed cutting, sawing and drilling of bone in the fields of Spine, ENT, Neuro and Maxillofacial Surgery.


ELAN[®] 4 Blades

Sagittal Saw Blades (Oscillating Saw Blades)




GA736

GA836



GA731

GA831





ELAN[®] 4 Blades


Sagittal Saw Blades (Oscillating Saw Blades)

5 sterile



GP491R	14 mm	13 mm	0.3 mm	0.3 mm	<p>» GA736, GA836</p>  
GP492R	15 mm	5 mm	0.3 mm	0.3 mm	
GP493R	20 mm	5 mm	0.3 mm	0.3 mm	
GP494R	20 mm	10 mm	0.3 mm	0.3 mm	
GP495R	20 mm	15 mm	0.3 mm	0.3 mm	
GP496R	25 mm	5 mm	0.3 mm	0.3 mm	
GP497R	25 mm	12 mm	0.3 mm	0.3 mm	


GP491R




GP496R



GP492R




GP497R




GP493R



GP494R



GP495R



Sagittal Saw Blades (Oscillating Saw Blades)




GP501R	15 mm	5 mm	0.5 mm	0.5 mm
GP502R	15 mm	10 mm	0.5 mm	0.5 mm
GP504R	20 mm	5 mm	0.5 mm	0.5 mm
GP505R	20 mm	10 mm	0.5 mm	0.5 mm
GP514R	25 mm	5 mm	0.5 mm	0.7 mm
GP515R	25 mm	10 mm	0.5 mm	0.7 mm
GP516R	25 mm	15 mm	0.5 mm	0.7 mm
GP518R	35 mm	10 mm	0.5 mm	0.5 mm
GP521R	35 mm	5 mm	0.5 mm	0.5 mm
GP522R	35 mm	10 mm	0.5 mm <td 0.7 mm	
GP523R	35 mm	15 mm	0.5 mm	0.7 mm
GP524R	35 mm	20 mm	0.5 mm	0.7 mm


» GA731, GA831




GP501R




GP502R




GP504R




GP505R



GP514R




GP515R




GP516R



GP518R




GP521R




GP522R



GP523R



GP524R



ELAN[®] 4 Blades

Sagittal Saw Blades (Oscillating Saw Blades)


1 sterile



GP535R	50 mm	10 mm	0.5 mm	0.8 mm	» GA731, GA831
GP536R	50 mm	15 mm	0.5 mm	0.8 mm	
GP537R	50 mm	20 mm	0.7 mm	0.9 mm	
GP538R	50 mm	25 mm	0.7 mm	0.9 mm	
GP539R	50 mm	30 mm	0.7 mm	0.9 mm	




GP535R




GP536R




GP537R



GP538R



GP539R



Sagittal Saw Blades (Oscillating Saw Blades)




GP511R	35 mm	10 mm	0.5 mm	0.5 mm	10 mm
GP526R	51 mm	15 mm	0.5 mm	6 mm	20 mm
GP527R	51 mm	15 mm	0.5 mm	7 mm	20 mm
GP528R	51 mm	15 mm	0.5 mm	8 mm	20 mm
GP529R	51 mm	15 mm	0.5 mm	9 mm	20 mm
GP530R	51 mm	15 mm	0.5 mm	10 mm	20 mm
GP531R	51 mm	15 mm	0.5 mm	11 mm	20 mm
GP532R	51 mm	15 mm	0.5 mm	12 mm	20 mm


» GA731, GA831



GP511R



GP526R -
GP532R




ELAN[®] 4 Blades

Reciprocating Saw Blades




GA737

GA837



GA732

GA832



Reciprocating Saw Blades




GP542R	13 mm	0.3 mm	0.3 mm
GP543R	20 mm	0.3 mm	0.3 mm
GP544R	13 mm	0.3 mm	0.3 mm
GP545R	20 mm	0.3 mm	0.3 mm

» GA737, GA837




GP542R



GP543R



GP544R





GP545R




ELAN[®] 4 Blades

Reciprocating Saw Blades




GP551R	84 mm	25 mm	0.40 mm	0.60 mm	» GA732, GA832  
GP561R	104 mm	25 mm	0.40 mm	0.60 mm	
GP552R	85 mm	27 mm	0.25 mm	0.35 mm	
GP553R	85 mm	27 mm	0.40 mm	0.60 mm	
GP562R	105 mm	27 mm	0.25 mm	0.35 mm	
GP563R	105 mm	27 mm	0.40 mm	0.60 mm	
GP554R	91 mm	33 mm	0.40 mm	0.60 mm	
GP564R	111 mm	33 mm	0.40 mm	0.60 mm	


GP551R



GP561R




GP552R
GP553R



GP562R
GP563R



GP554R



GP564R



Reciprocating Saw Blades



GP565R	115 mm	37 mm	0.4 mm	0.6 mm	
GP550R	73 mm	15 mm	0.4 mm	0.6 mm	
GP575R	79 mm	17 mm	0.4 mm	0.6 mm	4 mm
GP576R	86 mm	24 mm	0.4 mm	0.6 mm	6 mm

» GA732, GA832




GC615R	23 mm	0.25 mm	0.35 mm
GC633R	23 mm	0.25 mm	0.35 mm




GP559R	70 mm
GP569R	90 mm

GP565R




GP550R




GP575R



GP576R




GC615R



GC633R



GP559R



GP569R



ELAN[®] 4 Blades

Reciprocating Saw Blades

1 sterile




GP557R	70 mm	10.0 mm
GP558R	70 mm	12.5 mm
GP567R	90 mm	10.0 mm


» GA732, GA832




GP557R



GP558R




GP567R




ELAN[®] 4 Blades

Transversal Saw Blades



GA833

GA733



ELAN[®] 4 Blades

Transversal Saw Blades



GP581SU	8 mm	4 mm	0.3 mm	0.4 mm
GP583SU	10 mm	10 mm	0.5 mm	0.7 mm
GP584SU	15 mm	12 mm	0.5 mm	0.7 mm
GP587SU	11 mm	7 mm	0.6 mm	0.7 mm
GP588SU	11 mm	12 mm	0.6 mm	0.7 mm
GP591SU	8 mm	25 mm	0.5 mm	0.7 mm
GP594SU	8 mm	12 mm	0.5 mm	0.7 mm
GP595SU	14 mm	12 mm	0.5 mm	0.7 mm
GP598SU	37 mm	25 mm	0.5 mm	0.7 mm

» GA833, GA733



GP581SU **S Transversal**
to be used with GP580R

GP583SU **M Transversal**
to be used with GP582R & GP586R

GP584SU **M Transversal**
to be used with GP582R & GP586R

GP587SU **M Transversal**
to be used with GP582R & GP586R

GP588SU **M Transversal**
to be used with GP582R & GP586R

GP591SU **L Transversal**
to be used with GP590R

GP594SU **L Transversal**
to be used with GP590R

GP595SU **L Transversal**
to be used with GP590R

GP598SU **L Transversal**
to be used with GP590R



To be used with

GP580R	65 mm		S Transversal Saw Blades
GP582R	65 mm		M Transversal Saw Blades
GP586R	95 mm		M Transversal Saw Blades
GP590R	30 mm		L Transversal Saw Blades

ELAN[®] 4 Burrs



GA761 - GA764



GA861 - GA864



GA765 - GA769



GA865 - GA869



GA749




GA849



GB945R



GB941R



GB942R



GB943R




GB947R



ELAN[®] 4 Burrs

ELAN 4 Burrs - 1 Ring



Ø	0.6 mm	0.8 mm	1.0 mm	1.0 mm	1.4 mm	1.4 mm	1.8 mm	1.8 mm	2.3 mm	2.3 mm	2.7 mm	
I	GP111R	GP112R	GP113R	GP114R	GP115R	GP116R	GP117R	GP118R	GP120R	GP121R	GP122R	Rosen

Ø	2.7 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm	
I	GP123R	GP124R	GP125R	GP126R	GP127R	GP128R	GP129R	Rosen


Ø	3.0 mm	4.0 mm	5.0 mm	6.0 mm		Ø	4.0 mm	5.0 mm	6.0 mm	7.0 mm	
Soft Cut 											
I	GP133R	GP134R	GP135R	GP136R	Rosen	I	GP184R	GP185R	GP186R	GP187R	Twin-Cut Burr

Ø	0.6 mm	0.8 mm	1.0 mm	1.0 mm	1.4 mm	1.4 mm	1.8 mm	1.8 mm	2.3 mm	2.3 mm	2.7 mm	
I	GP141R	GP142R	GP143R	GP144R	GP145R	GP146R	GP147R	GP148R	GP149R	GP150R	GP151R	Diamond Burr

Ø	2.7 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm	
I	GP152R	GP153R	GP154R	GP155R	GP156R	GP157R	GP158R	Diamond Burr



ELAN® 4 Burrs - 1 Ring



Ø	2.0 mm	2.0 mm	2.3 mm	2.3 mm	3.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm	
	GP161R	GP162R	GP163R	GP164R	GP165R	GP166R	GP168R	GP169R	GP170R	Diamond Burr, Coarse

Ø	3.0 mm	4.0 mm	5.0 mm	6.0 mm	7.0 mm	
	GP173R	GP174R	GP175R	GP176R	GP177R	Diamond Burr, Extra Coarse

Ø	1.5 mm	2.0 mm	2.5 mm	3.0 mm	3.0 mm	
US Type 						
	GP201R	GP202R	GP203R	GP204R	GP205R	
	Neuro Cutter					



Ø	2.0 mm	2.5 mm	3.0 mm	
Aesculap Type 				
	GP208R	GP209R	GP210R	
	Neuro Cutter			

Ø	1.5 mm	1.5 mm	2.0 mm	2.0 mm	2.5 mm	2.5 mm	3.0 mm	3.0 mm	4.0 mm	4.0 mm	
	GP211R	GP212R	GP213R	GP214R	GP215R	GP216R	GP217R	GP218R	GP219R	GP220R	Neuro Cutter Diamond




ELAN[®] 4 Burrs

ELAN 4 Burrs - 1 Ring




1 sterile

Barrel Burr




1 sterile

Barrel Burr




1 sterile

Cone Burr




1 sterile

Acorn Burr




1 sterile

Oval Burr




1 sterile

Reverse Taper Burr Coarse Diamond




1 sterile

Rosen




1 sterile

Neuro Cutter




1 sterile

Lindemann




1 sterile

Twist Drill



1 sterile

Pin Cutter



ELAN® 4 Burrs - 2 Ring

Ø	2.3 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm	7.0 mm	8.0 mm	
II	GP301R	GP302R	GP303R	GP304R	GP305R	GP306R	GP307R	GP308R	GP309R	Rosen

Ø	2.3 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm	
II	GP311R	GP312R	GP313R	GP314R	GP315R	GP316R	GP317R	Diamond Burr

Ø	2.3 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm	
II	GP321R	GP322R	GP323R	GP324R	GP325R	Diamond Burr, Coarse

Ø	4.0 mm	5.0 mm	6.0 mm	
II	GP328R	GP329R	GP330R	Diamond Burr, Extra Coarse

Ø	2.0 mm	2.5 mm	3.0 mm	3.0 mm
US Type 				
II	GP361R	GP362R	GP363R	GP364R
	Neuro Cutter			



Ø	3.0 mm
Aesculap Type 	
II	GP367R
	Neuro Cutter

Ø	3.0 mm	3.0 mm
II	GP370R	GP371R
	Neuro Cutter Diamond	




ELAN[®] 4 Burrs

ELAN 4 Burrs - 2 Ring




1 sterile

Barrel Burr




1 sterile

Barrel Burr




1 sterile

Cone Burr




1 sterile

Acorn Burr




1 sterile

Oval Burr




1 sterile

Reverse Taper Burr Coarse Diamond




1 sterile

Craniotome Cutter - Spiral Type




1 sterile

Craniotome Cutter - Straight Type




1 sterile

Lindemann




1 sterile

Twist Drill



1 sterile

Twist Drill



ELAN® 4 Burrs – 2 Ring

Ø	1.0 mm	Ø	5.0 mm	6.0 mm	7.0 mm	8.0 mm	9.0 mm
II	GP350R	II	GP335R	GP336R	GP337R	GP338R	GP339R
	sterile		sterile				
Pin Cutter		Twin-Cut Burr					

Tungsten Carbide Burrs

Ø	2.3 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm	7.0 mm	8.0 mm
II	GP301TC	GP302TC	GP304TC	GP306TC	GP307TC	GP308TC	GP309TC
	sterile						
Rosen							

Ø	1.2 mm	1.6 mm	Ø	3.0 mm
II	GP395TC	GP396TC	II	GP367TC
	sterile			sterile
Taper Burr			Neuro Cutter	


Burrs for metal dissection

Ø	2.3 mm
II	GP399TC-SU
	sterile
Side Cutter	



ELAN[®] 4 Burrs


ELAN 4 MIS



GB920R



GB921R



10 cm



GB922R



GB925R



GB926R



13 cm



GB927R



ELAN® 4 MIS Burrs

Ø	3.0 mm	4.0 mm	Ø	3.0 mm	4.0 mm	Ø	2.0 mm	3.0 mm	Ø	2.0 mm	3.0 mm
L10	GP422SU	GP423SU	L13	GP442SU	GP443SU	L10	GP420SU	GP421SU	L13	GP440SU	GP441SU
	sterile			sterile			sterile			sterile	
Rosen Burr			Rosen Burr			Neuro Cutter			Neuro Cutter		

Ø	2.0 mm	3.0 mm	4.0 mm	Ø	2.0 mm	3.0 mm	4.0 mm	5.0 mm
L10	GP424SU	GP425SU	GP426SU	L13	GP444SU	GP445SU	GP446SU	GP447SU
	sterile				sterile			
Neuro Cutter Diamond				Neuro Cutter Diamond				

Ø	3.0 mm	4.0 mm	5.0 mm	Ø	3.0 mm	4.0 mm	5.0 mm	Ø	3.0 mm	4.0 mm	5.0 mm
L10	GP432SU	GP433SU	GP434SU	L13	GP452SU	GP453SU	GP454SU	L10	GP435SU	GP436SU	GP437SU
	sterile				sterile				sterile		
Diamond Burr				Diamond Burr				Diamond Burr Extra Coarse			

Ø	3.0 mm	4.0 mm	5.0 mm	Ø	3.0 mm	4.0 mm	5.0 mm	Ø	3.0 mm	4.0 mm	5.0 mm
L13	GP455SU	GP456SU	GP457SU	L10	GP428SU	GP429SU	GP430SU	L13	GP448SU	GP449SU	GP450SU
	sterile				sterile				sterile		
Diamond Burr Extra Coarse				Neuro Cutter Diamond Extra Coarse				Neuro Cutter Diamond Extra Coarse			



All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogs and price lists and to take legal actions.

©2019 AESCULAP. ALL RIGHTS RESERVED. PRINTED IN THE USA.
Aesculap is an equal opportunity employer

Aesculap, Inc. | 3773 Corporate Parkway | Center Valley, PA | 18034
Phone 800-282-9000 | Fax 610-791-6886 | www.aesculapusa.com

Aesculap, Inc. - a B. Braun company

DOC1360 Rev E 8/19