

New Extended Length Laparoscopic Instruments



Developed for challenging bariatric surgeries, and more.

Aesculap Laparoscopy

Part Number	New Instrument	Tip Description	Jaw Insert	Insulated Outer Tube	Standard Handle	Optional Handles
PO239R	Dolphin Nose Fixation Forceps 5 mm 420 mm	Fine tip attribute facilitates precise dissection and suture grasping	PO792R	PM981R	PO958R	PO959R PM950R PM954R PO960R PM951R PM957R PO961R PM953R
PO240R	Duck Bill Dissecting Forceps 5 mm 420 mm	Provides strong hold with a broad grip design and fine teeth	PO793R	PM981R	PO958R	PO959R PM950R PM954R PO960R PM951R PM957R PO961R PM953R
PO241R	Right-Angle Dissector 5 mm 420 mm	Features a serrated jaw for right-angle double action	PO794R	PM981R	PO958R	PO959R PM950R PM954R PO960R PM951R PM957R PO961R PM953R
PO242R	Overholt Dissector 90 Degree 10 mm 420 mm	Offers 90 degree curved jaws; suited for blunt dissection	PO795R	PM988R	PO958R	PO959R PM950R PM954R PO960R PM951R PM957R PO961R PM953R
PO192R	Universal Grasper 5 mm 420 mm	Fenestrated tip provides enhanced tissue gripping	PO796R	PM981R	PO959R	PO958R PM950R PM954R PO960R PM951R PM957R PO961R PM953R
PO245R	Atraumatic DeBakey Clamp Jaw 5 mm 420 mm - 40 mm Jaw	Features coarsely ribbed jaw for atraumatic tissue dissection	PO798R	PM981R	PO959R	PO960R PM951R PM957R PO961R PM953R PM958R PM950R PM954R

Aesculap, Inc. | 3773 Corporate Parkway | Center Valley, PA | 18034
Phone 800-282-9000 | Fax 610-791-6886 | www.aesculapusa.com

All rights reserved. Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogs and price lists and to take legal actions.

Aesculap, Inc. – a B. Braun company

©2020 AESCULAP. ALL RIGHTS RESERVED. PRINTED IN THE USA.
Aesculap is an equal opportunity employer

DOC1696 6/20