

# Ligation Clip System

A Reliable and Cost-Effective System

OPERATE WITH GREATER PRECISION



AESCULAP®

# Ligation Clip System

A Reliable and Cost Efficient System that Delivers Superior Performance



## Titanium Clips

- Positive grip atraumatic design
- Consistent closing without scissoring
- Chevron shape provides a tip-to-tip closure encompassing the vessel and tissue



## Appliers for Open Surgery

- Available in 6", 8" and 11" lengths; 20°, 25° and 65° angle
- Color-coded to match popular clip sizes
- Consistent applier springs minimize inadvertent squeezing of jaws
- Spring unlatches for easy cleaning
- Fine-tip models available



## Appliers for Laparoscopic Surgery


- Available in 30 cm, 33 cm, 45 cm
- Rotation knob allows for one finger, 360° rotation of applier shape
- Modular design allows for disassembly during cleaning

## Ligation Clips Cartridges

- Easy identification through color-coding
- Pull-tab adhesive feature on cartridge bottom sticks where you want it
- Six sizes available:



J1180-1  
Small  
(yellow)



R1180-1  
Small-Wide  
(red)




L5180-1  
Small-Medium  
(purple)



B2180-1  
Medium  
(blue)



V3120-1  
Med.-Large  
(green)



O4120-1  
Large  
(orange)


Ordering Information						
Titanium Ligation Clips	Small	Small-Wide	Small-Medium	Medium	Medium-Large	Large
Box of 6-clip Cartridges	<b>J1180-1</b> (30 cartridges/box)	<b>R1180-1</b> (30 cartridges/box)	<b>L5180-1</b> (30 cartridges/box)	<b>B2180-1</b> (30 cartridges/box)	<b>V3120-1</b> (20 cartridges/box)	<b>O4120-1</b> (20 cartridges/box)
Box of 24-clip Cartridges	<b>J1120-1</b> (5 cartridges/box)	<b>R1120-1</b> (5 cartridges/box)		<b>B2120-1</b> (5 cartridges/box)		
Open Appliers						
6" (15 cm) 25° angle	<b>PJ115-1B</b>	<b>PR115-1B</b>		<b>PB215-1B</b>		
8" (20 cm) 20° angle	<b>PJT120-1B</b>	<b>PRT120-1B</b>	<b>PLT520-1B*</b>	<b>PBT220-1B*</b>		
8" (20 cm) 25° angle					<b>PV320-E</b>	<b>PO420-1B</b>
8" (20 cm) 65° angle					<b>PVC320-E</b>	<b>POC420-E</b>
11" (28 cm) 25° angle				<b>PB228-1B</b>	<b>PV328-E</b>	<b>PO428-1B</b>
11" (28 cm) 65° angle				<b>PBC228-1B</b>	<b>PVC328-E</b>	<b>POC428-1B</b>
Laparoscopic Appliers						
30 cm	<b>LCASLS1R</b> (5.5 mm dia)		<b>LCASLS2R</b> (5.5 mm dia)	<b>LCASLS3R</b> (6.5 mm dia)		
33 cm				<b>LCASLS4R</b> (10 mm dia)	<b>LCASLS5R</b> (10 mm dia)	<b>LCASLS6R</b> (12.5 mm dia)
45 cm						<b>LCASLS5-45R</b> (10 mm dia)

**Indications for Use:**

Hemostatic clip is designed for the intended use of ligating blood vessels. The clip has been specially designed to insure occlusion of vessels and prevent any slippage once applied. The clip has applications in many types of surgical procedures where hemostasis is required or radiographic marking is necessary.

Choose the size of clip to fit the procedure making certain the tissue to be occluded fits completely within the clip. The Vitalitec International, Inc., clip can be left in vivo without sequela as it is biologically inert.

\* Denotes T-Style Open Applier



Through collaborative  
excellence we will improve the  
quality of a patient's life and  
meet the needs of the changing  
healthcare environment.

Technical alterations are possible. This leaflet may be used for no other purposes than offering, buying and selling of our products. No part may be copied or reproduced in any form. In the case of misuse we retain the rights to recall our catalogs and price lists and to take legal actions.

©2021 AESCULAP, INC. ALL RIGHTS RESERVED. PRINTED IN THE USA.  
Aesculap is an equal opportunity employer

Aesculap, Inc. | 3773 Corporate Parkway | Center Valley, PA | 18034  
Phone 800-282-9000 | Fax 610-791-6886 | [www.aesculapusa.com](http://www.aesculapusa.com)

Aesculap, Inc. - a B. Braun company

DOC426 Rev G 1M 12/21

